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Aluminum Salt Stain Kit (Lillie Method)

Cat: G3020

Size: 3×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent		3×50mL	Storage
Reagent (A): Aluminum Stain Solution	A1: Lillie Buffer	50mL	2-8°C
	A2: Aluminum Reagent	1g	RT, avoid light
Before use, mix A1 with A2 as the ratio of 10mL:0.2g to form Aluminum Stain Solution. It is ready to use.			
Reagent (B): Lillie Decolorizing Solution		50mL	2-8°C
Reagent (C): Fast Green Solution		50mL	RT, avoid light

Introduction

Most of aluminum in the body is excreted with feces through the digestive tract, and a small part is accumulated in testicles, kidneys, spleen, muscles, bones and brain tissues. The accumulation of aluminum in the brain can cause the degeneration of brain nerves, memory decline, intelligence and personality will also be affected, and even present Alzheimer's disease. When the amount of aluminum in the body is more than 5-16 times of the normal amount, it can inhibit the absorption of phosphorus in the intestine and interfere with the normal metabolism of calcium and phosphorus in the body.

The principle of Aluminum Salt Stain Kit (Lillie Method) is that aluminum reagent combines with aluminum ion to form red complex.

Protocol (for reference only)

1. Conventional fixation, wax immersion and embedding. Dewax to distilled water.
2. Preparation of Aluminum Stain Solution: Before use, mix A1 with A2 as the ratio of 10mL:0.2g to form Aluminum Stain Solution. It is ready to use and can be stored at 4 °C for a short time.
3. Preheat the Aluminum Stain Solution to 80-85 °C in advance, pour it into the plastic slide dyeing container, immerse the slide into the Aluminum Stain Solution and place it in the microwave oven, and heat it for 30s in the second level (120 °C). Take out the dyeing container and put it at room temperature for 10min. Wash with distilled water for 3 times.
4. Place in Lillie Decolorizing Solution for 5s. Wash with distilled water for 3 times.
5. Re-dyeing with Fast Green Solution for 3 min. Wash with distilled water for 3 times.
6. Erect and air dry. Slightly immerse in xylene and seal with resinene.

Result

Aluminum	Red
Background	Green

Note

1. Pay attention to the power of the microwave oven when heating. If the power is too large, the heating time shall be reduced.
2. For your safety and health, please wear experimental clothes and disposable gloves.