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Cat: G3535

Size: 5×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

| Reagent | Size | Storage |
|---|------|--------------------|
| Reagent(A): Gomori Oxidant | 50mL | RT, avoid light |
| Reagent(B): Oxalic Acid Solution | 50mL | RT |
| Reagent(C): Ferric Ammonium Sulfate Solution | 50mL | RT, avoid light |
| Reagent(D): Gomori Ammoniacal Silver Solution | 50mL | 2-8°C, avoid light |
| Reagent(E): Gomori Reductant | 50mL | RT |

Introduction

Reticulin Stain Kit(Gomori Sliver Method)mainly goes through the steps of oxidation, bleaching, mordant dyeing, silver dipping, reduction, and re-dyeing.The difference between this method and the modified Gordon-Sweets method is that Gordon-Sweets method uses acid oxidant and nuclear fast red re-dyeing solution. Frozen section, low temperature section and collodion section can be used for reticular fiber staining. All kinds of fixatives can be used, and the heavy metal mercury salt or osmium salt fixatives occasionally produce some non-specific silver background. It is often used to distinguish the nature and source of tumor, cancer and sarcoma, lymphosarcoma and reticulosarcoma, hemangioendothelioma and hemangiopericytoma, etc.

Protocols(for reference only)

1. Fix the tissue in 10% formalin fixative, then conventionally dehydrate and embed.
2. Cut into paraffin section in 4μm thick, and dewax to distilled water.
3. Put the section on the dyeing frame, drop in Gomori Oxidant and oxidize for 5min.
4. Wash with water slightly.
5. Bleach in Oxalic Acid Solution for 1-2min.
6. Rinse in running water for 2min and slightly with distilled water.
7. Mordant dyeing with Ferric Ammonium Sulfate Solution for 5min.
8. Wash slightly with tap water and distilled water.
9. Add Gomori Ammoniacal Silver Solution and dye for 3min.
10. Wash slightly with distilled water.
11. Reduce by Gomori Reductant for 1min.
12. Rinse in running water for 10min.
13. Conventionally dehydrate and transparent, seal with resinene.

Result

| | |
|-----------------|---------------------------|
| Reticular Fiber | Black |
| Collagen Fiber | Yellow or Yellowish Brown |
| Nucleus | Brown or Black Brown |

Note

1. Glass container must be soaked in washing solution for one day, washed with tap water and washed twice with distilled water.
2. 10% formalin fixative is a suitable fixative. It is not suitable to use fixative containing mercury such as Zenker's Solution, otherwise it is easy to cause nonspecific deposition of sections.
3. Gomori Ammoniacal Silver Solution is not very stable and sensitive to light. It should be kept away from light at 4 °C and used after restoring to room temperature.
4. For your safety and health, please wear experimental clothes and disposable gloves.