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Movat-Russell Modified Pentachrome Stain Kit

Cat: G3701

Size: 10×50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	10×50mL	Storage	
Reagent(A): Hypo Solution	50mL	RT	
Reagent(B): Alcian Blue Staining Solution	50mL	2-8°C, avoid light	
Reagent(C): Alkaline Alcohol Solution	50mL	RT	
Reagent(D): Weigert Hematoxylin Solution	D1: WH Souldution A	30mL	RT, avoid light
	D2: WH Souldution B	20mL	RT, avoid light
	D3: WH Souldution C	10mL	RT, avoid light
Before use, mix D1, D2 and D3 as the radio of 3:2:1 to form Weigert Hematoxylin Solution. It is not suitable to prepare in advance.			
Reagent(E): WH Differentiation	50mL	RT, avoid light	
Reagent(F): Biebrich Scarlet – Acid Fuchsin Solution	E1: BS-AF Solution A	40mL	RT, avoid light
	E2: BS-AF Solution B	10mL	RT, avoid light
Before use, mix E1 with E2 as the radio of 4:1 to form Biebrich Scarlet-Acid Fuchsin Solution. It is not suitable to prepare in advance.			
Reagent(G): Phosphotungstic Acid Solution	50mL	RT, avoid light	
Reagent(H): Weak Acid Differentiation Solution	50mL	RT	
Reagent(I): Alcoholic Safran Solution	50mL	RT, avoid light	
Reagent(J): Russell Mordant Solution	50mL	RT, avoid light	

Introduction

In a narrow sense, connective tissue refers to three kinds of fibers: collagen fiber, reticular fiber and elastic fiber. There are many staining methods for connective tissue, such as Masson Trichrome Staining, Van Gieson Staining, Gomori Ammoniacal Silver Staining, Mallory Phosphotungstate Hematoxylin Staining. However, the above staining methods only focus on one or two kinds of tissue staining. Movat-Russell Modified Pentachrome Staining is popular for its rich and bright colors. This staining method is mainly used to show atherosclerotic plaque. Weigert Hematoxylin Solution is used to dye the nucleus, Biebrich Scarlet-Acid Fuchsin Solution is used to dye the cytoplasm, Alcoholic Safran Solution is used to dye collagen tissue, Alcian Blue Staining Solution dye matrix (proteoglycan). Because of the complexity of the operation of the kit, its dyeing effect has a lot to do with the operator's experience and quantity, so it is not easy to dye the kit with satisfactory results at the same time.

Self Provided Materials

Series of ethanol, Distilled water, Microwave oven

Protocols(for reference only)

1. For paraffin section, conventionally dewax and dehydrate in series of ethanol.
2. Take a proper amount of Russell Mordant Solution and put it into microwave oven, heat it for 30-60s moderately, and treat the section immediately for 10mins. Rinse in running water for 10mins. (See Note 3)
3. Treat the section in Hypo Solution for 5min and rinse in distilled water 2-3 times for 30s each.
4. Stain the section in with Alcian Blue Staining Solution for 20mins and rinse in distilled water for 30-60s.
5. Preheat the Alkaline Alcohol Solution in a water bath or oven at 45-60 °C, and then treat the section in the Alkaline Alcohol Solution for 10min. Rinse in distilled water for 2-5min.
6. Add the section in prepared Weigert Hematoxylin Solution and dye avoid light for 10-30min.
7. Slightly wash in running water and rinse in distilled water 2-3 times for 30s each.
8. Decolor the section with WH Differentiation for 10s, rinse with distilled water for 10s.
9. Add the section in prepared Biebrich Scarlet-Acid Fuchsin Solution and dye avoid light for 1min. Rinse in distilled water 2-3 times for 30s each.
10. Treat the section in Phosphotungstic Acid Solution for 3-5min and directly treat in Weak Acid Differentiation

Solution for 3min. Rinse in distilled water 2-3 times for 30s each.

11. Dehydrate: 95% ethanol for 1min, 100% ethanol 2 times for 1min each .
12. Stain with Alcoholic Safran Solution for 3min.
13. Dehydrate: absolute ethanol 2 times each for 1min. Transparent in xylene and seal.

Result

Nucleus and Elastic Fiber	Deep Purple to Black
Collagen Protein and Reticular Fiber	Red
Proteoglycan	Blue
Cellulose	Dark Red
Cardiac Muscle and Smooth Muscle	Red

Note

1. Because of staining power and tissue section, all five colors may not be displayed after staining. Pay attention to anti stripping treatment.
2. The thickness of section is about 5 μ m, and the differentiation of elastic fiber is usually completed within 2-3min.
3. It is very important to wash with water to remove the Weigert Hematoxylin Solution. If the washing is not enough, the subsequent dyeing steps will be inhibited.
4. This method can show cryptococcus neoformans and dye them bright blue.