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Cat:G1274

Size:3×50mL

Storage:RT, avoid light, valid for at least 1 year.

Kit Components

Reagent	3×50mL	Storage
Reagent A: Ziehl-Neelsen Fuchsin Solution	50mL	RT , avoid light
Reagent B: Ziehl-Neelsen Destaining Solution	50mL	RT
Reagent C: Methylene Blue Solution	50mL	RT , avoid light

Introduction

The cell walls of mycobacteria contain a large amount of lipid in the outside of the peptidoglycan, so it is not easy staining. When the mycobacterium acid of mycobacteria combines with dye, it is difficult to be destained by acidic destaining solution. This method is called Acid-Fast Bacillus Stain. Traditional dyeing methods need to promote staining by heating and extending the time. The most representative method is Ziehl-Neelsen method, which is recommended in the WHO and China Tuberculosis Control Program.

Acid-Fast Bacillus(AFB)Stain Kit(Ziehl-Neelsen's Method) is a kind of thermal dyeing solution. The dyeing principle is that under heating condition, the mycobacterium acid combines with Fuchsin Stain and binds to a complex. Mycobacterium remains red after Methylene Blue re-dyeing, and the others are blue. The kit is more suitable when the effect of cold dyeing is not good.

Protocol(for reference only)

1. Pick up the samples to be examined with the inoculation ring and coat them on the glass slide. Then heat and fix the section.
2. Add Ziehl-Neelsen Fuchsin Solution to the section and heat it slightly with flame until vapor appears. Generally, the dyeing process should be at least 5 min (if necessary, add Ziehl-Neelsen Fuchsin Solution to prevent the evaporation of dyeing solution). Rinse with distilled water.
3. Destain with Ziehl-Neelsen Destaining Solution until there is no red, usually for 1 min. Rinse with distilled water.
4. Dyeing with Methylene Blue Solution for 1 min. Rinse with distilled water.
5. Gently absorb water and naturally dry. View under the oil immersion microscope.

Result

Acid-fast bacteria	Red
Background	Blue

Note

1. Cover the reagent bottle tightly after each use to prevent reagent evaporation and contamination.
2. If it is difficult to achieve flame heating in operation step 2. A small amount of fuchsin dye solution can be separated in advance into the EP tube, and preheated at 80-100 °C for 5-10min after sealing. Add 80-100 °C hot water to the incubation box for thermal insulation dyeing (add dye solution if necessary to prevent the dye solution from drying).
3. If the staining solution is used for paraffin section, the solvent such as gasoline turpentine must be used to replace xylene for transparency, dewaxing and other steps to avoid the integrity of the lipid layer on the surface of the ring-breaking bacteria and reduce the positive rate.
4. The above reagents are irritating to human body. Please pay attention to appropriate protection.
5. For your safety and health, please wear experimental clothes and disposable gloves.