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Cat :G1580

Size:100T

Storage:-20°C,avoid light, valid for 1 year.

Kit Components

| Reagent | 100T | Storage |
|--|-------|-------------------|
| Reagent(A): β -Gal Fixative | 100mL | 2-8°C,avoid light |
| Reagent(B): X-Gal Solution | 5mL | -20°C,avoid light |
| Reagent(C): β -Gal Stain Solution A | 1mL | 2-8°C,avoid light |
| Reagent(D): β -Gal Stain Solution B | 1mL | 2-8°C,avoid light |
| Reagent(E): β -Gal Stain Solution C | 100mL | 2-8°C |
| Mix Reagent B,C,D,E, in proportion of 5:1:1:93 as Dyeing Working Solution. It is ready to use. | | |

Introduction

The Senescence β -Galactosidase Staining Kit is a staining kit for senile cells or tissues based on the up-regulation of SA- β -Gal(senescence-associated β -galactosidase) activity during aging. Aging of cells or tissues can be observed under optical microscopes. Cell aging is also considered as a way for organisms to inhibit tumors and a potential cause of organism aging. The kit can be used for senescence detection of cultured cells and tissue sections.

Using X-Gal as substrate, the Senescence β -Galactosidase Staining Kit produces dark blue products catalyzed by aging-specific β -galactosidase. Cells or tissues expressing β -galactosidase that turn blue are easily observed under optical microscopy. This kit only stains senescent cells. For tissue sections or blocks, the number of samples that can be detected depends on the size of the sample. For ordinary sections or 6-well plate, at least 100 samples can be detected.

Protocol(for reference only)

Note: It is recommended to use polypropylene (PP) container to prepare dyeing working solution, but polystyrene (PS) container cannot be used to prepare dyeing working solution.

The incubation at 37 °C during dyeing could not be carried out in carbon dioxide incubator.

Adherent Cells Stain

1. For cells cultured in 6-well plates,remove the cell culture medium and wash once with PBS or HBSS. Add 1 ml β -Gal Fixative to each pore and fix at room temperature for 15 min. For other types of culture plates, the amount of fixative solution and subsequent solution is operated according to this ratio.
2. Remove β -Gal Fixative and wash cells 3 times with 1×PBS for 3 min each time. (*priority step see note 4*)
3. Remove washing solution and add 1 ml Dyeing Working Solution to each hole.
4. When incubate overnight at 37°C, the 6-well plate can be sealed with parafilm to prevent evaporation.
5. View under the optical microscope. If you can't view and count in time, you can remove the Dyeing Working Solution, add 2ml PBS, and store it at 2-8 °C for several days; moreover, it can be stored for a long time after sealing with aqueous gelatin sealing agent.

Suspension Cell Stain

1. Collect the cells into 1.5ml eppendorf tube by centrifugation, wash once with PBS or HBSS, and fix at room temperature for 15 min with 1 ml β -Gal Fixative. Slowly shake on the shaking table to avoid forming cell clumps.
2. Centrifuge and remove β -Gal Fixative, wash cells with 1×PBS 3 times for 3 min each time. (*priority step see note 4*)
3. Centrifuge and remove the washing solution then add 0.5-1 ml Dyeing Working Solution to each tube. Incubate overnight at 37°C.
4. Drip some stained cells onto section and view under optical microscope. If not view in time, centrifuge and remove the dyeing working solution then adding 1 ml PBS at 2-8°C can keep for several days. If centrifuged, cells can also be taken for smears, and sealed with water based sealing agent, it can be stored for a long time at 2-8°C.

Tissue Sections Stain

1. For frozen sections, follow the following steps directly.
2. Add a suitable volume of β -Gal Fixative to cover the tissue adequately and fix it for more than 15 min. (*priority step see note 4*)
3. Wash the tissues 3 times with PBS for more than 5 min each. Then incubate with washing solution for 3 min.
4. Remove washing solution and add appropriate amount of Dyeing Working Solution.
5. When incubate overnight at 37°C, it can be placed in the wet-box to prevent evaporation.
6. View under optical microscope. If not view in time, can seal with glycerol gelatin for storage.

Note

1. The β -Gal Fixative has certain corrosiveness and toxicity. Please pay attention to the protection during operation.
2. The dyeing reaction of β -Gal depends on the specific pH condition and incubation should not be in a carbon dioxide incubator. The higher concentration of carbon dioxide in the incubator used for cell culture will affect the pH value of the dyeing working solution and lead to dyeing failure.
3. The precipitation in the reagent B will not affect staining. The precipitation would dissolve completely after full mix or vortex. As a routine, the precipitation must be dissolved and mixed before use.
4. The priority step is using reagent E: β -Gal Stain Solution C to wash and incubate instead of 1 \times PBS before staining.
5. PP containers or glass containers should be used in the preparation of dyeing working solution instead of PS containers. But dyeing can be carried out in PS containers, such as 6-well plate.
6. If long-term storage is needed, the excess liquid can be poured out and dried, and then sealed with neutral gum.
7. For your safety and health, please wear experimental clothes and disposable gloves.