Davidson 固定液

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规格: 100mL/500mL

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Storage: RT, valid for 1 year.

Introduction

Davidson's Fixative is a fast fixative, especially suitable for the preparation of tumors, bone marrow samples, gynecological materials, fatty breast and tissue biopsy. It can also be used for testicular biopsy and overnight fixation of fatty lymph nodes. In addition, The optimum choice of Davidson's Fixative is the whole eye fixative, which can provide the best tissue preservation conditions and avoid retinal detachment.

Davidson's Fixative does not contain any dyes (such as eosin), mercury or other metals.

Protocol (for reference only)

1. Fix the larger sample or several smaller samples with a fixative of about 100ml. Take out the sample and remove it into the fixative immediately.
2. Fix the sample in fixative for 1-2 days. The fixation of small samples should not exceed 24 h.
3. After fixation, the sample can be stored in formalin or 70% ethanol.

Note

1. During fixation, the tissue should be completely submerged to prevent uneven fixation.
2. The fixing time changes due to different thickness of tissue samples. If it is not used for a long time (interval of more than 1 month). It is recommended to store it in a cool and dark place.
3. If the specimen is especially hard (such as eyes), it can be stored in 70% ethanol for one day after fixation, then transferred to 90% ethanol for one day, and finally stored in absolute ethanol for one day, so as to make the specimen easier to cut.
4. For your safety and health, please operate carefully in ventilated environment and wear experimental clothes and disposable gloves.