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## Methenamine Silver Plating Stain Kit For Basement Membranes

Cat: G1790

Size: 6×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent		6×50mL	Storage
Reagent(A): Oxidant		50mL	2-8°C, avoid light
Reagent(B): Ammonium Ferric Sulfate Solution		50mL	RT, avoid light
Reagent(C): Methenamine Silver Working Solution	C1: Methenamine Silver Solution	25mL	2-8°C, avoid light
	C2: Sodium Borate Solution	25mL	RT
Before use, mix C1 with C2 in equal amount to form Hexamine Silver Working Solution. It is ready to use.			
Reagent(D): Gold Chloride Solution		50mL	2-8°C, avoid light
Reagent(E): Hypo Solution		50mL	RT
Reagent(F): Light Green Solution		50mL	RT, avoid light

### Introduction

Methenamine Silver Staining is a classical method to display basement membrane. After oxidation, the mucopolysaccharide in the basement membrane was exposed to aldehyde group, which reduced methenamine silver to black metallic silver. Gold Chloride Solution can transform silver into a more stable metal gold and make the background clearer. Methenamine Silver Staining can clearly show the basement membrane, which is widely used in nephrosis. It is often used to observe the morphological changes of glomerular capillary basement membrane in inflammatory injury, such as rupture, proliferation, folding and so on.

### Protocol (for reference only)

1. Dewax the paraffin section in xylene and rehydrate in ethanol.
2. Add Oxidant onto the section and oxidize for 15mins. Wash with distilled water for 2min.
3. Stain with Ammonium Ferric Sulfate Solution for 10mins. Wash with distilled water for 2min.
4. Mix C1 and C2 in equal amounts to prepare Methenamine Silver Working Solution, which must be used within 2 hours.
5. After thoroughly cleaning the slices, pour out excess water, add newly prepared Methenamine Silver Working Solution dropwise to cover the tissue, and dye at 60 °C water bath or constant temperature for 20-30 min until the slices turn tobacco yellow or black. Wash in distilled water for 1 min.
6. Soak in Hypo Solution for 1 min. Wash with distilled water for 1mins.
7. Color with Gold Chloride Solution for 1-2mins. Wash with distilled water for 2min.
8. Re-dyeing with Light Green Solution for 1 min. Wash with distilled water for 1min.
9. Dehydrate in 95% ethanol, transparent by xylene and seal with resinene.

### Result

Basement membrane of glomerular capsule and glomerular capillaries	Reddish Brown to Black
Background	Green

### Note

1. The glass container used in the experiment shall be soaked in the washing liquid in advance, rinsed and dried.
2. The prepared silver hexamine working solution is disposable and cannot be put for a long time. It is recommended to use it within 2 hours.
3. When gold chloride is used for color matching, it is necessary to observe and control the color matching results under the microscope.
4. The light green solution is a lining reagent, which can be replaced with other lining reagents according to the observation needs, such as G1100- Eosin Y Stain Solution, For HE, G1120-Hematoxylin-Eosin (HE) Stain Kit, etc.
5. For your safety and health, please wear experimental clothes and disposable gloves.