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## Modified Sirius Red Stain Kit (No Picric Acid)

**Cat:** G1472

**Size:** 2×50mL/2×500mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent		2×50mL	2×500mL	Storage
Reagent(A):Iron Hemat-oxylin Working solution	A1:Iron Hematoxylin Stock Solution	25mL	250mL	RT, avoid light
	A2:Iron Hematoxylin Diluent Solution	25mL	250mL	RT, avoid light
Before use, mix A1 with A2 in equal to form Iron Hematoxylin Working Solution, it is recommended to use within 4 h.				
Reagent(B):Sirius Red Staining Solution		50 mL	500 mL	RT, avoid light

### Introduction

Collagen fiber is the most widely distributed and abundant fiber in connective tissue. It is widely distributed in various organs. Type I collagen fibers are mainly bone, skin and tendon; Type II collagen fibers are mainly cartilage collagen; Type III collagen fibers are mainly found in embryonic tissues, adult blood vessels and gastrointestinal tract; Type IV collagen fibers are mainly in the basement membrane. Sirius red and its lining dye solution are strong acid dyes, which are easy to combine with the basic groups in collagen molecules and have strong adsorption. Polarized light microscopy shows that the collagen fibers have the property of positive uniaxial birefringence. When combining with Sirius red compound staining solution, it could enhance the birefringence and improve the resolution, so as to distinguish the two types of collagen fibers.

Modified Sirius Red Stain Kit (No Picric Acid) is mainly composed of Iron Hematoxylin Stain and Sirius Red Stain, without picric acid. It is mainly used in the study of collagen fiber abnormality or fiber proliferation in various tissue diseases. Under the ordinary optical microscope, the collagen fibers of heart and blood vessels and other tissues are dyed red, which is helpful for the classification and grading of various fibrotic diseases under the polarized light microscope. Immunohistochemical technique can also display type I and III collagen fibers, but the antibodies used are expensive and time-consuming, while Sirius red staining reagent is cheap and easy to operate.

### Self Provided Materials

10% formalin fixative solution, ordinary optical microscope or polarized light microscope

### Protocol(for reference only)

1. Fix the tissue in 10% formalin fixative solution and embed by routine dehydration.
2. Cut the slices into 3-8 μm. Conventional dewaxing to water.
3. Prepare Iron Hematoxylin Staining Solution before use, stain the section with it for 5-10min, and wash it with distilled water for 10-20s to remove the excess staining solution.
4. (Optional) Immerse in tap water for 5-10min, and wash with distilled water for 3 times, 5-10s each time.
5. Sirius Red Staining Solution shall be used for dripping dyeing for 15-30min. For tissues that are easy to be stained, the dyeing time can be controlled within 5-10min.
6. Rinse slightly with running water to remove the staining liquid on the surface of the slice.
7. Dehydrate by series of ethanol and transparent by xylene, finally seal with neutral gum.

### Result

1) Optical microscope observations	Collagen fiber	Red
	Nucleus	Brown to Black
	Muscle fiber	Yellow
2) Polarized light microscope observations	Type I collagen fiber	Strong Orange Yellow or Light Red
	Type III collagen fiber	Green

### Note

1. In order to display clearly under the polarized light microscope, the slice thickness should be 6-7μm.
2. The nucleus without iron hematoxylin staining may be light red, which is a normal phenomenon, but it will not occur after iron hematoxylin staining.
3. This method is stable and not easy to fade.
4. For your safety and health, please wear experimental clothes and disposable gloves.