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Plant Tissue Stain Kit (PAS-NYS method)

Cat: G4809

Size: 3×50mL

Storage: 2-8°C, avoid light, valid for 6 month.

Kit Components

Reagent	3×50mL	Storage
Reagent (A): PAS Oxidant	50mL	2-8°C, avoid light
Reagent (B): Schiff Reagent	50mL	2-8°C, avoid light
Reagent (C): NYS Solution	50mL	RT, avoid light

Introduction

Starch in plant tissues and sugars in cell walls are exposed to aldehyde groups under the action of oxidants. The aldehyde groups combine with Schiff reagent to produce a red reaction, which is called PAS staining. Proteins in plants are often stored in the form of aleurone particles in plant cells, in a solid state, surrounded by a layer of membrane, forming spherical particles. Naphthol yellow S forms a yellow complex with proteins. Therefore, PAS staining combined with naphthol yellow can simultaneously display starch granules and proteins in plant tissues.

This product is used for staining protein and starch granules in plant tissues. After staining, the plant starch granules are purple red, the cell walls are purple red, and the proteins are yellow.

Protocol(for reference only)

1. Take fresh plant tissue for routine fixation, recommend 10% formalin, and perform routine dehydration and embedding.
2. Conventional dewaxing of paraffin slices to distilled water; Freeze the slices directly into distilled water and soak for 2min.
3. Stain with PAS Oxidant for 10-15min, wash twice with distilled water for 5-10s each time.
4. Stain with Schiff Reagent in dark at room temperature for 10-20min, wash twice with distilled water for 5-10s each time.
5. Stain with NYS Solution for 5-10min, then rinse quickly with distilled water for 5-10s.
6. Dehydrate by series of ethanol and each for 3-5s, transparent with xylene and seal with neutral gum.

Result

Starch granules	Purple red granules
Cell wall	Purple red
Protein	Yellow

Note

1. Slice dewaxing should be as clean as possible, otherwise it will affect the staining effect.
2. The PAS Oxidant can be reused, and the dye solution is transparent. If it shows obvious yellow color, impurities, or tissue staining is too light, a new dye solution needs to be replaced.
3. Schiff Reagent needs to be stored at 2-8°C, and before use, it should be taken out and restored to room temperature. When the color of the solution turns noticeably red, it is recommended to replace it with a new dye solution.
4. For your safety and health, please wear laboratory clothes and disposable gloves for operation.