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Tissue Transparent Dewaxing Solution, Benzene Free

Cat: G2150

Size: 100mL

Storage:RT, valid for 3 years.Used up within 6 months after opening.

Introduction

This product is a colorless and transparent liquid with a slight natural plant fragrance, which can be used for micro sectioning of plant, animal, and human anatomy. This product is basically insoluble in water and has moderate volatility. Long term immersion in tissue will not make it brittle, and it can also enhance the wax immersion effect, which is beneficial for subsequent slicing.

It is mainly composed of oils extracted from natural plants and saturated alkane derivatives, and does not contain volatile toxic substances such as xylene. It is suitable for the treatment of pathological and cytological samples, and has good compatibility with alcohol and paraffin. It is an ideal substitute for xylene in biological tissue transparency treatment, and has excellent transparency effect.

Protocol (for reference only)

1. Transparency of the organization: Transfer the dehydrated tissue block into a transparent solution, and depending on the size of the tissue block, it should be transparent for 1-4 hours. To achieve complete transparency, it is recommended to replace the dewaxed transparent solution once in the middle.
2. Dewaxing of slices: Preheated paraffin slices are immersed in transparent liquid, and depending on the thickness of the slices, dewaxing is performed for 10 to 30 minutes. To fully dewaxing, it is recommended to replace the dewaxing transparent liquid at least once in the middle.
3. Transparent before sealing: After dehydration with ethanol in the slicing series, enter the transparent solution twice, and let it dry slightly after 1 minute each time to seal the film. This product is a non-polar solvent and may not be compatible with some neutral gum products when sealing. Before use, it is recommended to take 20ul (about 1 drop) of sealing agent and mix it with an equal amount of this product for observation. If it is cloudy or has obvious precipitation, it cannot be directly used for transparency before sealing. (See Note 2)

Note

1. When the tissue becomes transparent after dehydration, the amount of transparent liquid should not be too small, which is 5-10 times the total volume of the tissue block. The tissue block should be thoroughly dehydrated before becoming transparent.
2. This product is mainly composed of saturated alkane derivatives, which may cause white wax precipitation when mixed with some market products containing xylene sealing gum. If it is cloudy, it can be washed with xylene before sealing or used with benzene free sealing gum, such as our company G8591/G8592/8593 -Benzene Free Neutral Gum (B-C-D type).
3. The biochemical reagents produced by our company are basically non sterile packaging without special labeling. If they are used for cell experiments, please prepare them in advance.
4. Product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
5. For your safety and health, please wear White coat and disposable gloves.