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规格: 1mL(100T)/5mL(500T)

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**Cat:**A7631

**Size:**1mL(100T)/5mL(500T)

**Storage:**2-8 °C, avoid light, valid for 1 year.

### Introduction

Alamar Blue Cell Viability Assay Reagent provides a simple, rapid, reliable and safe method for cell proliferation and cytotoxicity detection, which is suitable for high-throughput test. The main component of the reagent is a redox indicator. In the oxidation state, it is purple blue without fluorescence, but in the reduction state, it turns into a reduction product with pink or red fluorescence. The absorption peak is 530-560nm, while the scattering peak is 590nm.

In the process of cell proliferation, the ratio of NADPH / NADP, FADH / FAD, FMNH / FMN and NADH / NAD increase and the cell is in a reducing environment. The dyes absorbed in the cells were reduced by these metabolic intermediates and cytochromes, then released to the cells and dissolved in the culture medium, which changed the culture medium from non fluorescent indigo blue to fluorescent pink. Finally, it can be detected by ordinary spectrophotometer or fluorescence photometer. The absorbance and fluorescence intensity are in direct proportion to the number of active cells.

### Protocol(for reference only)

1. Add 10  $\mu$ l of detection reagent into 100  $\mu$ l cell suspension, incubate in cell incubator for 2-6 h, and enter the next step if the color of culture medium change from indigo blue to pink.
2. It is recommended to use a fluorescent enzyme reader for detection. Make the excitation wavelength between 530-560 nm, the emission wavelength as 590 nm, and record the relative fluorescence unit (RFU).
3. Draw standard curve or cell growth curve: vertical coordinate (Y axis) is relative fluorescence unit (RFU); horizontal coordinate (X axis) is cell number or time point or drug concentration.

### Note

1. Appropriate density of cells can increase the sensitivity of detection. For 96 well plate, we suggest to inoculate 100  $\mu$ l cells per well which own the cell concentration range for adherent cells in 100-10000 / well, suspension cells in 2000-50000 / well, and the culture medium is blank control.
2. For 384 well plates, the cell concentration and inoculation amount are halved.
3. The whole process should be aseptic, because microbial pollutants can also reduce the detection reagent and affect the experimental results.