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规格: 50mL/100mL/500mL

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[1] Dan Xiao, Yue Zhang, Rui Wang, et al. Emodin alleviates cardiac fibrosis by suppressing activation of cardi-ac fibroblasts via upregulating metastasis associated protein 3. Acta Pharmaceutica Sinica B. April 2019. (IF 5.808)

[2] Chunping Yang, Tianxing Lv, Bin Wang, et al. The Damaging Effects of Pedunsaponin A on Pomacea canaliculata Hemocytes. Toxins. July 2019. (IF 3.895)

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Storage: 2-8°C , valid for 1 year.

Introduction

Trypan Blue, is a cell reactive dye, which is often used to detect the integrity of cell membrane and the survival rate of cells. It is one of the most commonly used dead cell identification methods in tissue and cell culture. Normal living cells have intact cell membrane structure and can reject trypan blue, so the cells will not be stained blue. However, cells with loss of activity or incomplete cell membrane can be stained blue by trypan blue. It is generally considered that the cell has died when the integrity of the cell membrane is lost. Trypan Blue Stain Solution, 0.4% is most commonly used in microscopy for cell counting.

Trypan blue can be engulfed by macrophages, so it can be used as a live stain for macrophages. Apoptotic bodies also have trypan blue rejection phenomenon. After Trypan blue staining, cell viability can be accurately quantified by counting directly under the microscope or after photographing under a microscope. Trypan blue stain solution, 0.4% is most commonly used for microscopic cell counting, which only takes 3-5 minutes and is easy to operate.

Protocol(for reference only)

Collecting Cells

For adherent cells, first digest them with trypsin; For suspended cells, can collect directly. Centrifuge the collected cells at 1000-2000rpm for 3 minutes, discard the supernatant, add an appropriate amount of PBS to resuspend the cells, and prepare a single-cell suspension.

Trypan Blue Stain

Mix 1 part of Trypan Blue Stain Solution, 0.4% and 9 part of cell suspension (cell dilution).

Incubate the mixture approximately for 3mins at room temperature.

Cell Count

Take a drop of the trypan blue / cell mixture to a hemacytometer. Place the hemacytometer on the stage of a binocular microscope and focus on the cells. Non-viable cells will be stained blue. Count the unstained (viable) and stained (nonviable) cells separately in the hemacytometer.

$$\text{viable cells (\%)} = \frac{\text{total number of viable cells permL of aliquot}}{\text{total number of cells permL of aliquot}} \times 100$$

Note

1. The staining time should not be too long, otherwise the living cells will gradually accumulate dye and stain, which will make the detection result low.
2. If there is precipitation in the dyeing solution before dyeing, it should be filtered out and reused.
3. There is a potential cancer risk. For your safety and health, please wear laboratory clothes and disposable gloves.