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## Trypan Blue Stain Solution, Ready To Use

**Cat:** C0041

**Size:** 100mL

**Storage:** 2-8°C , valid for 1 year.

### Introduction

Trypan Blue Stain Solution, 0.4% is an azo dye, derived from toluidine. It is also known as diamine blue and Niagara blue. It is widely used as a vital stain to distinguish the viable cells from the non-viable cells. It selectively stains the dead tissues or cells blue, while live cells or tissues with intact cell membrane are not colored and are very selective in the compounds that pass through the membrane. Among them, apoptotic bodies also have anti staining properties, while macrophages can ingest trypan blue through endocytosis.

This product is a low concentration ready to use trypan blue staining solution, using an improved buffer system with strong stability and less precipitation. By counting directly under the microscope or taking photos under the microscope, the cell survival rate can be accurately quantified. The dyeing time is only 3-5 minutes, and the operation is simple.

### Protocol(for reference only)

#### Collecting Cells

For adherent cells, first digest them with trypsin; For suspended cells, can collect directly. Centrifuge the collected cells at 1000-2000rpm for 3 minutes, discard the supernatant, add an appropriate amount of PBS to resuspend the cells, and prepare a single-cell suspension.

#### Trypan Blue Stain

Mix 9 part of Trypan Blue Stain Solution, Ready To Use and 1 part of cell suspension (cell dilution).

Incubate the mixture approximately for 3min at room temperature.

#### Cell Count

Take a drop of the trypan blue / cell mixture to a hemacytometer. Place the hemacytometer on the stage of a binocular microscope and focus on the cells. Non-viable cells will be stained blue. Count the unstained (viable) and stained (nonviable) cells separately in the hemacytometer.

$$\text{viable cells (\%)} = \frac{\text{total number of viable cells permL of aliquot}}{\text{total number of cells permL of aliquot}} \times 100$$

### Note

1. The staining time should not be too long, otherwise the living cells will gradually accumulate dye and stain, which will make the detection result low.
2. If there is precipitation in the dyeing solution before dyeing, it should be filtered out and reused.
3. There is a potential cancer risk. For your safety and health, please wear laboratory clothes and disposable gloves.