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Cat: G0170

Size: 100T

Storage: -20°C, avoid light, valid for 1 year.

### Kit Components

Reagent	100T	Storage
Reagent A : MDC Stain Solution	1mL	-20°C, avoid light
Reagent B : 10×Wash buffer	20mL	2-8°C
Reagent C : Collection buffer	10mL	2-8°C

### Introduction

Autophagy is a process that cells are stimulated to phagocytize their own cytoplasm or organelles and eventually degrade the phagocytes in lysosomes. Autophagosome is a circular or elliptic structure encapsulated by bilayer membranes. It contains cytoplasm, longevity proteins, complex of abnormal proteins, damages or redundant organelles, viruses and bacteria.

Dansylcadaverine (MDC) is a fluorescent pigment and an eosinophilic stain. It is usually used to detect the formation of autophages as a specific marker dye. Its detection excitation filter wavelength is 355 nm and the blocking filter wavelength is 512 nm. The Autophagy/Cytotoxicity Dual Stain Kit(MDC Method) is suitable for autophagy staining of cultured cells and can counterstain with EB.

### Self Provided Materials

Fluorescence microscope, Low-speed centrifuge, EB stain

### Protocol(for reference only)

#### For Suspension Cell Stain

1. Dilute 10× Wash buffer to 1× with distilled water.
2. 800g centrifugation for 5 min. Collect the cell in the tube and wash once with 1×Wash buffer of 300-400 μL. Remove the supernatant.
3. Add a certain amount of 1×Wash buffer to suspend cells, count and regulate the cell concentration to 10<sup>6</sup>/mL.
4. Take 90μL cell suspension into a new EP tube and add 10 μL MDC stain solution and gently mix well. Incubate for 15-45 min at room temperature and avoid light. (See Note 1)
5. 800g centrifugation for 5 min. Collect the cell in the tube and wash twice with 1×Wash buffer of 300-400 μL. Remove the supernatant.
6. Add 100μL Collection buffer to suspend the cell, drop it onto the slide and cover the slide.
7. View under the fluorescence microscope (blocking filter wavelength is 512 nm), count and photograph.

#### For Adherent Cell Stain

1. remove the culture medium and wash with Diluted 1× Wash buffer.
2. Dilute MDC Stain Solution with 1× Wash buffer by 1:10, then add to cultivation plate and stain for 15-45min.
3. Wash with 1×Wash buffer twice and 3min for each time.
4. View under the fluorescence microscope (blocking filter wavelength is 512nm), count and photograph.

### Result

Normal cells	Cells are evenly dyed yellow-green fluorescence.
Apoptotic cells	Collapsed chromatin and punctate nucleus are dyed into dense green particles of varying sizes.

### Note

1. The combination of MDC stain solution and EB reagent can distinguish normal cells, apoptotic cells, and necrotic cells. It has certain toxicity, please handle with caution. If combined staining is required, a final concentration of 0.2 can be added to the cell suspension with MDC staining solution μ M's EB staining solution, followed by normal operation.
2. Here, 10 times diluted MDC staining solution can be added according to the principle of reducing the volume of the culture medium added to the well plate by half. The dilution ratio can also be adjusted according to the number of cells to avoid false positives caused by high concentration.
3. During the operation, pay attention to reducing the exposure time of reagents to strong light.
4. For your safety and health, please wear experimental clothes and disposable gloves.