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规格：4×10mL/4×100mL/4×500mL

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Cat: G1060

Size: 4×10mL/4×100mL/4×500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	4×10mL	4×100mL	4×500mL	Storage
Reagent A: Gram's Crystal Violet Solution	10mL	100mL	500mL	2-8°C
Reagent B: Gram's Safranin Solution	10mL	100mL	500mL	2-8°C
Reagent C: Gram's Iodine Solution	10mL	100mL	500mL	2-8°C, avoid light
Reagent D: Gram's Decolorizing Solution	10mL	100mL	500mL	RT

Introduction

The Gram Stain is a different staining technique most widely applied in microbiology. Gram staining is based on the ability of bacteria cell wall to retaining the crystal violet dye during solvent treatment. The cell walls for Gram-positive microorganisms have a higher peptidoglycan and lower lipid content than gram-negative bacteria. Bacteria cell walls are stained by the Gram's Crystal Violet Solution. Iodine is subsequently added as a mordant to form the crystal violet-iodine complex so that the dye cannot be removed easily. However, subsequent treatment with Gram's Decolorizing Solution dissolves the lipid layer from the gram-negative cells. As a result, the diffusion of the violet-iodine complex is blocked, and the bacteria remain stained. The length of the decolorization is critical in differentiating the gram-positive bacteria from the gram-negative bacteria. Finally, counterstain with Gram's Safranin Solution, the gram-negative bacteria is stained a pink color.

Protocol (for reference only)

1. Prepare a bacterial smear on clear, dry glass slide, and fix it through a gentle flame 3-4 times while moving the slide in a circular fashion to avoid localized overheating. It is appropriate not to burn your hands.
2. Gram Staining:
 - 1) Flood with Gram's Crystal Violet Solution for 1 min. Wash with distilled water.
 - 2) Flood the smear with Gram's Iodine Solution for 1 min. Wash with distilled water.
 - 3) Decolorize with Gram's Decolorizing Solution for 20 to 60 s until the blue dye no longer flows.
 - 4) Wash with distilled water. Counterstain with Gram's Safranin Solution for 1 min. Wash with distilled water.
 - 5) Allow the slide to air dry or blot dry between sheets of clean bibulous paper and view under oil immersion lens.

Result

Gram-Positive Organisms	Bluish Purple
Gram-Negative Organisms	Pinkish Red

Note

1. The specimen smear should not be too thick and should be carried out in strict accordance with the operation requirements. If the smear is thick, the decolorization time should be extended until purple no longer appears.
2. The temperature of slide on flame should not be too high. If the iodine solution becomes colorless, change it.
3. If there is no oil microscope observation condition or the slices need to be preserved for a long time, the neutral resin seal can be added after drying and observed under 40 times microscope.
4. When washing, the action should be gentle, and the washing bottle should be used along the diagonal direction of the slide to avoid washing off the bacteria.