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## Crystal Violet Ammonium Oxalate Solution, 0.4%

**Cat:** G1070

**Size:** 100mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Introduction

Crystal Violet Ammonium Oxalate Solution, 0.4% is a component of Gram's staining method. It is mainly composed of ammonium oxalate, crystal violet and ethanol. In Gram-positive bacteria staining, cells are stained with ammonium oxalate crystal violet and Gram iodine solution to form an insoluble complex. The complex could not pass through the cell wall and could not be destained easily, so it remained purple. After dyeing, the bacteria formed a sharp contrast with the environment, and the morphology, arrangement and some structural characteristics of bacteria could be clearly observed for classification and identification.

### Self provided materials

Inoculation ring or other tools for picking up bacteria, Alcohol lamps, Slides, Optical microscope.

### Protocol (for reference only)

Operate according to the specific requirements of the experiment or refer to the following methods of Gram staining.

1. **Smear:** Take the bacterium to be examined and smear it in the center of the slide or add a little sterile water on the slide. Mix the bacterium with water evenly and smear it in a thin layer.
2. **Drying:** After smear, it can be dried naturally at room temperature or heated slightly on alcohol lamp to make it dry quickly.
3. **Fixation:** Hold the slide with sample facing upward and move 3-5 times quickly outside the flame of alcohol lamp for 1 second each time. The temperature should not be too high to prevent the degeneration of bacterial proteins. Stain after cooling. It can also be fixed with methanol or ethanol.
4. **Primary dyeing:** Add Crystal Violet Ammonium Oxalate Solution, 0.4% to dye for 1-2 min, then rinse the dye solution with water.
5. **Mordant dyeing:** Add Gram iodine solution, and cover the slide at room temperature for 1-2 min, then wash it with water.
6. **Destaining:** Add the destaining solution and shake it for 10-30 s until the destaining solution does not appear purple. Rinse the destaining solution with water immediately to stop the reaction.
7. **Re-dyeing:** Re-dyeing with Safranin stain solution for 30-60 s, then washing with water.
8. **Observation:** View under the Oil immersion microscope after drying.

### Result

Gram-positive bacteria	Blue to Purple
Gram-negative bacteria	Red

### Note

1. Before smear, make the circle mark on the back in advance to judge the position of subsequent test.
2. When taking bacteria, pay attention to self-protection. When pulling or plugging the test tube plug, we should burn the test tube orifice slightly through the flame, and finally burn the inoculation ring on the flame to sterilize.
3. When heating the fixed smear, pay attention to the slide not too close to the flame. Generally, the temperature of the slide should not over 60°C.
4. The key of Gram staining is to strictly control the destaining degree, and the destaining time should be judged by experience. If destaining is excessive, positive bacteria can be mistakenly stained as negative bacteria; if destaining is insufficient, negative bacteria can be mistakenly stained as positive bacteria.
5. The culture time of the bacteria to be examined will affect the staining. When the culture time of the positive bacteria is too long or the bacteria has died or dissolved, the result is often negative.
6. For your safety and health, please wear experimental clothes and disposable gloves.