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# Mayer's Hematoxylin Stain Solution, For IHC

**Cat:** G1080

**Size:** 10mL/100mL/500mL

**Storage:** RT, valid for at least 1 year.

## Introduction

Hematoxylin, as an alkaline dye and one of the two main staining reagents for HE staining, is usually used for nuclear staining. It can be roughly divided into two types according to different dyeing operations and applicable samples. One is the progressive dyeing hematoxylin with less hematoxylin content, which is suitable for lining dyeing and usually does not need differentiation. The other is the degenerative staining hematoxylin with high content of hematoxylin and clear and long-lasting staining but needs differentiation.

This product is a progressive hematoxylin staining solution, which is suitable for counterstaining of tissue sections in immunohistochemistry. Usually, it does not need differentiation, and the nuclei are blue after staining. This product is working fluid and can be used directly. The dyeing solution can be reused for many times, but the dyeing effect will gradually decrease.

## Protocol(for reference only)

1. Incubate in Mayer's Hematoxylin Stain Solution (For IHC) for 30s-2min or 5-10min depending on desired signal intensity. Wash with distilled water for 1min to remove the excess staining solution.
2. (Optional) Differentiate with differentiation solution(0.05% -0.5% hydrochloric acid solution) for a few seconds, if the color is too deep. The staining results were observed under a washing microscope with rapid distilled water.
3. Wash with tap water for 10min or other bluing solution for 3-10min to fully blue.
4. According to the type of color rendering solution, seal the tablet directly with water sealing agent. Or dehydrate in gradient alcohol, transparent by xylene and seal with resinene.

## Result

Nucleus	Blue
---------	------

## Note

1. When using this kit for the first time, it is recommended to take 1-2 samples for pre-test.
2. Staining, differentiation, and blue-turning time all require pre-experimental optimization to be optimal, and the differentiation steps are optional.
3. Preserved at room temperature and protected from light, and the best temperature is 2-8 °C, which is at least one year.
4. Light staining is recommended in the process of staining. Generally, it is only necessary to be able to distinguish the nucleus. Too deep color may affect the color of the cytoplasm.
5. For your safety and health, please wear experimental clothes and disposable gloves.

## Related Products

G1100 Eosin Y Stain Solution, For HE

G1120 Hematoxylin-Eosin(HE) Stain Kit

G1865 Scott's Bluing Solution

P1031 PBST, 1× (pH7.2-7.4, 0.01M)

P1110 Paraformaldehyde, 4%

S2100 Mounting Medium, antifading

S2130 Mounting Medium, Antifading(PVP)

C1010 Sodium Citrate buffer(0.01mol/L, pH6.0, powder)