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规格: 2×50mL/2×100mL/2×250mL

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**Storage:** RT, valid for 1 year.

### Kit Components

Reagent	2×50mL	2×100mL	2×250mL	Storage
Reagent A: Crystal Violet Solution	50mL	100mL	250mL	RT, avoid light
Reagent B: Copper Sulfate Solution	50mL	100mL	250mL	RT

### Introduction

Most bacterial capsules are composed of polysaccharide. The capsule differs from the slime layer that most bacterial cells produce in that it is a thick, detectable, discrete layer outside the cell wall. Bacterial capsules are non-ionic, so neither acidic nor basic stains will adhere to their surfaces. Therefore, the best way to visualize them is negative stain.

For examination under the microscope, stain the bacteria and background darker than the capsule, which doesn't stain. Capsular material is very moist (slimy) and any heating will cause it to shrink. So do not heat-fix. Also, heating may cause the bacterial cell to shrink resulting in a clear zone around the cell - which may cause cells which don't have capsules to appear.

### Protocol(for reference only)

1. Prepare a smear of *Capsulate bacteris*, allow the smear to dry and then fix. (**Note : DO NOT heat or blot dry!**)
2. Drip an appropriate amount of Crystal Violet Solution onto the section (to cover the bacterial membrane), stain for 5-7 min, and drain excess staining solution to air dry.
3. Add an equal amount of Copper Sulfate Solution, shake slightly to differentiate, and then pour it away. Do not rinse with running water. Rinse moderately to ensure a purple background. Blot dry with bibulous paper. then add 1-2 drops cedar wood oil to prevent crystal formation in ice cream.
4. View under oil immersion lens.

### Result

Background	Blue to Purple
Thalli	Purple
Capsul	Colorless or Pale Purple

### Note

1. The degree of staining can be adjusted appropriately by changing the staining time, preferably with similar bacterial bodies and backgrounds.
2. During decolorization, do not rinse with water. Reagent (B) must be used to avoid excessive background decolorization affecting observation.
3. The section shall not be heated, dried or fixed to avoid washing the glass with water and preventing the capsule from shrinking or losing.
4. Please wear the experimental clothes and disposable gloves.