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规格：3×50mL/3×100mL/3×250mL

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Cat: G1170

Size: 3×50mL/3×100mL/3×250mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×50mL	3×100mL	3×250mL	Storage
Reagent A: Acid-Fast Stain Solution A	50mL	100mL	250mL	RT, avoid light
Reagent B: Acid-Fast Stain Solution B	50mL	100mL	250mL	RT
Reagent C: Acid-Fast Stain Solution C	50mL	100mL	250mL	RT, avoid light

Introduction

Acid Fast Stain is for the demonstration of Mycobacterium tuberculosis in tissue sections and smears. It is often used for preliminary identification of Mycobacterium. Compared with ordinary bacteria, the cell wall of Mycobacterium contains more lipid, so it is difficult to dye and must be heated to be stained by Carbol-Fuchsin Solution. Mycobacterium can resist the decolorization of acid decolorizing solution such as hydrochloric acid alcohol, while other bacteria and cells are decolorized. After re-dyeing with basic methylene blue, the Mycobacterium is still red, and other bacteria and background are blue.

Protocol(for reference only)

1. Perform a bacterial smear, and fix by flame after drying.
2. Drop Acid-Fast Stain Solution A on the smear and cover the sample. Heat the smear on the alcohol lamp, do not boil, remove it in case of steam, and add dye if necessary to avoid drying up. Heat for 3-5 min.
3. Move away fire and stand for 5 min, then wash in distilled water.
4. Differentiate in Acid-Fast Stain Solution B until no color runs off the smear. Complete decolorization avoids producing false-positive results.
5. Wash in distilled water for 3-5 min. Re-dyeing with Acid-Fast Stain Solution C for 30s.
6. Rinse in tap water and air dry. View the smear under oil immersion lens.

Result

Acid-Fast Bacilli	Red
Non-Acid-Fast Bacilli	Blue
Background	Blue

Note

1. The dyeing degree can be adjusted by changing the dyeing time or temperature.
2. After each reagent is used up, please cover it quickly to avoid volatilization.
3. This product has certain irritancy and corrosiveness. Please wear experimental clothes and disposable gloves when using.