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
## Modified Oil Red O Stain Kit, Environmental Friendly

**Cat:** G1263

**Size:** 4×50mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	4×50mL	Storage
Reagent(A):Oil Red O Buffer	50mL	RT
Reagent(B):Oil Red O Staining Solution	50mL	RT, avoid light
Reagent(C):Oil Red O Differentiation Solution	50mL	RT
Reagent(D): Mayer Hematoxylin Solution	50mL	RT, avoid light

### Introduction

Lipid is the general name of neutral fat, lipids and their derivatives. Its common physical characteristics are insoluble in water and soluble in organic solvents (such as ethanol, acetone, etc.). Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and other methods are often used for staining of lipids in tissues. The traditional method uses sudan dyes. Recently, it has been found that azo dye oil red O is more suitable for fat dyeing. Oil red O is a strong fat dye, which is easy to combine with neutral fat and has a slightly poor binding force with phospholipids. Its dyeing principle is generally considered that fat is dyed by similar miscibility or physical adsorption. The solubility of dyes in intracellular lipids is greater than that in original solvents. Therefore, when dyeing, dyes are transferred from organic solvents to lipids to dye fats.

Modified Oil Red O Stain Kit(Polyol Method) uses an odorless and low volatile environment-friendly solvent to replace classical solvent with certain risks, which greatly reduces the use risk and operation difficulty without affecting the location and size of lipids, and is very friendly to experimental operators. The result of fat positive staining is orange red to red, and the specific color depends on the type and concentration of lipids.

### Protocol(for reference only)

1. After sampling, fix tissues with G2161-Neutral buffered Formalin,10% or P1110-Paraformaldehyde, 4% precooled by 2-8 °C refrigerator for full fixation or quick freezing with liquid nitrogen to prepare fresh frozen sections. Generally, the section is recommended to be cut in 10-12um.
2. Before staining, rewarm the slices to room temperature. Drop the Oil Red O Buffer onto the section for 2-5min. (see Note 1)
3. Pour out the excess buffer without washing, and directly add the preheated Oil Red O Staining Solution to cover the section and dye at room temperature for 15min.
4. Pour out the excess staining solution and add the same amount of Oil Red O Differentiation Solution for differentiation treatment for 2-5min until the section color is uniform.
5. Soak in distilled water for 1min to fully remove the differentiation solution, and stain with Mayer Hematoxylin Solution for 2min.
6. After washing with distilled water for 1min, blue with tap water for 5-10min until the nucleus is clear blue.
7. Use the preheated and melted S2150-Glycerol Gelatin aqueous slide mounting medium to seal the section and observe it under the optical microscope.

### Result

Neutral Fat	Orange or Red
Phospholipid	Pink Red
Nucleus	Blue

### Note

1. Lipids are easily soluble in organic solvents. Generally, frozen sections are used for staining. It is recommended to use precooled G2161-Neutral buffered Formalin,10% or P1110-Paraformaldehyde,4% for fixation to prepare frozen sections. At the same time, it is not recommended to use acetone for secondary fixation of quick-frozen tissue sections.
2. Reagent(A):Oil Red O Buffer has the functions of fixation and buffering. For quick-frozen sections, the incubation time (5min) can be appropriately prolonged to properly fix the tissues. If the tissue is sectioned after fixation, the incubation time can be appropriately shortened (2min).

3. The dye solution and buffer solution are viscous. It is recommended to suck them after cutting the mouth of the suction head, or preheat them properly before use in winter (25-28 °C). At the same time, it is normal that it cannot be completely dumped because it is viscous. You can directly follow the operating steps for the next step.
4. The staining results cannot be stored for a long time, so they should be observed and photographed as soon as possible.
5. For your safety and health, please wear lab clothes and disposable gloves.