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## Acid-Fast Bacillus(AFB)Stain Kit(Modified Kinyoun's Method)

Cat:G1273

Size:3×50mL

Storage:RT, avoid light, valid for 1 year.

### Kit Components

Reagent	3×50mL	Storage
Reagent(A): Modified Kinyoun Fuchsin Solution	50mL	RT, avoid light
Reagent(B): Kinyoun Destaining Solution	50mL	RT
Reagent(C): Methylene Blue Solution	50mL	RT, avoid light

### Introduction

The cell walls of mycobacteria contain a large amount of lipid in the outside of the peptidoglycan, so it is not easy staining. When the mycobacterium acid of mycobacteria combines with dye, it is difficult to be destained by acidic destaining solution. This method is called Acid-Fast Bacillus Stain. Traditional dyeing methods need to promote staining by heating and extending the time. The most representative method is Ziehl-Neelsen method, which is recommended in the WHO and China Tuberculosis Control Program.

Acid-Fast Bacillus (AFB) Stain Kit (Modified Kinyoun's Method) belongs to cold dyeing solution, which needs no heating and is relatively safe for acid-fast dyeing solution. The dyeing principle is that at room temperature, Mycobacterium acid and carbolic acid fuchsin combine to form a complex. After methylene blue dyeing, Mycobacterium is still red, while others are blue. This method is lighter than Kinyoun's method and more suitable for weak acid fast dyeing.

### Self Provided Materials

Inoculation Ring, Slide, Distilled Water, Microscope

### Protocol(for reference only)

1. Pick up the samples to be examined with the inoculation ring and coat them on the glass slide. Then dry it in air.
2. Drop Modified Kinyoun Fuchsin Solution to stain for 5 min. Rinse with distilled water without heating.
3. Destain with Kinyoun Destaining Solution until there is no red for 3min. Rinse with distilled water.
4. Dyeing with Methylene Blue Solution for 3-5min. Rinse with distilled water.
5. Gently absorb water and naturally dry.
6. View under the oil immersion microscope.

### Result

Acid-Fast Bacteria or Weak Acid Fast Bacteria	Red
The Others	Light Blue

### Note

1. Cover the reagent bottle tightly after each use to prevent reagent evaporation and contamination.
2. If the staining solution is used for paraffin section, the solvent such as gasoline turpentine must be used to replace xylene for transparency, dewaxing and other steps to avoid the integrity of the lipid layer on the surface of the ring-breaking bacteria and reduce the positive rate.
3. The above reagents are irritating to human body. Please pay attention to appropriate protection.
4. For your safety and health, please wear experimental clothes and disposable gloves.