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Cat: G1316

Size: 100mL/500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

Neutral red acts as a pH indicator, changing from red to yellow between pH 6.8 and 8.0. The molecular formula is $C_{15}H_{16}N_4 \cdot HCl$, and the molecular weight is 288.78. In neutral or slightly alkaline environment, living cells of plants can absorb a large number of neutral red and excrete into vacuoles.

Neutral Red Solution, 0.1% can be used as a red nuclear counterstain in various histological procedures.

Protocol (for reference only)

1. Sample Treatment

For paraffin section, dewax in xylene twice for 5-10mins each. Rehydrate with series ethanol (100%, 90%, 80%, 70%) for each level 3mins, and finally in distilled water 3min.

For frozen section restore to RT and wash with distilled water for 2mins.

For cultured cell, fix with 4% PFA for more than 10mins. Then wash with distilled water twice for 2mins each.

2. Neutral Red Solution Staining

- For the above treated samples, stain with Neutral Red Solution, 0.1% for 0.5-1h.
- Rinse briefly in distilled water or distilled water.
- View under the microscope and photograph.

Note

- If dehydration, transparency and sealing treatment are needed, xylene, Series Ethanol and resinene or other sealing agent shall be prepared by oneself.
- When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test.
- When staining with Neutral Red Solution, the time can be adjusted according to the dyeing results and requirements.
- For your safety and health, please wear experimental clothes and disposable gloves.