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Cat: G1320

Size: 100mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Nuclear Fast Red Solution is commonly used in tissue section staining, and the nucleus is red after staining. This product is working solution and can be used directly. The dye solution can be reused many times, but the dyeing effect will gradually decrease.

Protocol(for reference only)

This product may produce suspended solids or a small amount of precipitation due to flocculation. It is recommended to take supernatant or cool it after boiling water bath for 5-10min. (see Note 2)

1. For paraffin section, dewax to water, and proceed other dyeing. After dyeing, wash with distilled water for 2-5mins.
2. Immerse the tissue section in Nuclear Fast Red Solution, 0.1%, or directly drop it on the sample to dye for 5-10mins (light staining nucleus). The dyeing time can be adjusted according to the dyeing degree.
3. Rinse in tap water or distilled water to remove excess dye.
4. Conventional dehydration and transparency, then seal with resinene.

Note

1. Section dewaxing should be as clean as possible.
2. This product is a colloidal solution, which is stored at low temperature (lower than 25 °C) or stored for a long time. Suspended solids or a small amount of precipitation are generated due to flocculation, which is a normal phenomenon and generally does not affect the use. If the colloid solution is evenly dispersed in the boiling bath, tighten the bottle cap for 5-10min to recover the turbid solution.
3. When using this reagent for the first time, it is recommended to take 1-2 samples for pretest.
4. Light staining is recommended in the dyeing process, which can usually distinguish the nucleus. Too deep color may affect the cell color.
5. For your safety and health, please wear experimental clothes and disposable gloves.