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**Cat:** G1355

**Size:** 3×50mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	3×50mL	Storage
Reagent A: Mallory Stain A	50mL	RT , avoid light
Reagent B: Mallory Stain B	50mL	RT , avoid light
Reagent C: Mallory Stain C	50mL	RT , avoid light

### Introduction

Trichrome Stain usually refers to the staining of nucleus and the selective stain of collagen and muscle fibers. The dyeing principle of this method is related to the size of anionic dye molecules and the permeability of the tissues: the size of the molecules is reflected by the molecular weight, and the small molecular weight can easily penetrate the compact and low permeability tissues, while the large molecular weight can only enter the loose and high permeability tissues.

Mallory Trichrome Stain Kit can dye collagen fiber blue, nucleus and muscle fiber red, and red blood cell orange.

### Protocol(for reference only)

1. Dewax paraffin sections and rehydrate in graded of alcohol.
2. Stain with Mallory Stain A for 3-5min.
3. Stain with Mallory Stain B for 1-2min.
4. Stain with Mallory Stain C for 5-10min.
5. Wash by 95% alcohol for 10-15s.
6. Conventional dehydration, transparentize by xylene.
7. Seal with resinene.

### Result

Red blood cell	Orange
Muscle fiber	Red
Collagen fiber	Blue

### Note

1. In the dyeing process, go directly to the next step after 2/3/4 steps without washing.
2. If there are too many slices to be operated in time, and the reagents can not be contaminated with each other at the same time, the weak acid solution or 0.5-2% glacial acetic acid aqueous solution can be used between 2/3/4 steps for cleaning and temporary storage. The staining results can also be observed with a small amount of weak acid under the microscope.
3. It is recommended to use potassium dichromate fixative or bouin fixative to fix tissues. If the tissue is easy to be over-stained, it can be sensitized with bouin fixative or G1349-trichromatic staining mordant before dyeing after the paraffin section is rehydrated.
4. For your safety and health, please wear experimental clothes and disposable gloves.