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## Saffron-O and Fast Green Stain Kit (For Plant)

**Cat:** G1375

**Size:** 2×50mL/2×100mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	2×50mL	2×100mL	Storage
Reagent(A): Saffron-O Staining Solution	50mL	100mL	RT, avoid light
Reagent(B): Fast Green Staining Solution	50mL	100mL	RT, avoid light

### Introduction

Saffron-O is an alkaline dye which can show the lignified, corked and keratinized tissues of vascular plants. Fast Green is an acid dye, which is widely used in dyeing cells and plant tissues. The vessels can be stained with Saffron-O and the sieve tubes can be stained with Fast Green.

Saffron-O and Fast Green Stain Kit (For Plant) is suitable for dyeing roots, stems and leaves of higher plants.

### Protocol(for reference only)

1. Fix the tissue in 4% paraformaldehyde solution and dehydrate routinely.
2. Dewax to water and then stain with Saffron-O Staining Solution for 1-2 h.
3. Wash with tap water to remove the excess dye.
4. Treat with 50%, 70% and 80% gradient alcohol for 1-2min.
5. Dye in Fast Green Staining Solution for 1 min.
6. Wash away excess dye with anhydrous ethanol
7. Transparent with xylene and seal with resinene.

### Result

Lignified cell wall and nucleus	Red
Parenchyma cell, cytoplasm	Green

### Note

1. For small tissue samples (stem with diameter less than 1.5 mm or leaf with thickness less than 1 mm), negative pressure treatment is not necessary.
2. Saffron-O and Fast Green Stain is stable and can be used repeatedly for a long time.
3. For your safety and health, please wear lab clothes and disposable gloves.