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Prussian Blue Iron Stain Kit (With Nuclear Fast Red)

Cat: G1422

Size: 2×50mL/2×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent		2×50mL	2×100mL	Storage
Reagent(A):Perls Stain Working Solution	A1:Perls Stain A	25mL	50mL	RT, avoid light
	A2:Perls Stain B	25mL	50mL	RT
Before use, mix equal parts of A1 and A2 to form Perls Stain. It is not suitable to prepare in advance.				
Reagent(B): Nuclear Fast Red Solution		50mL	100mL	RT, avoid light

Introduction

The Prussian Blue Reaction is a method of staining non heme iron (containing heme, iron transporters, iron deposition, etc.) in cells or tissues containing granular iron. It plays an important role in the validation of diseases such as ischemic anemia, hemolytic anemia, hemochromatosis, and abnormal hemoglobin metabolism. This staining method can effectively distinguish between hemosiderin and other pigments, but due to the principle of staining, it is only suitable for staining tissues and cells with a high content of trivalent iron ions. In combination with other enhancing reagents, it can also be used for the detection of iron death models.

The Prussian blue reaction was first proposed by Perls and is therefore most commonly referred to as the Perl iron staining method. However, since its proposal in 1867, numerous chemists and pathologists have optimized this reaction, among which the Lillie method, which is more widely used and well optimized, can label divalent and trivalent iron separately. This kit is a ready to use Lillie method ferric iron staining kit, prepared with high-quality raw materials, with stable reagents and good staining specificity.

Self Provided Materials

Fixative(See Note 2), Series of ethanol, Distilled water

Protocol(for reference only)

1. The organization recommends fixation in neutral formalin fixative (10%) , conventional dehydration and embedding. (See Note 2)
2. The recommended thickness for paraffin sections is 5µm. Conventional dewaxing to water. The recommended thickness for frozen sections is 10µm. Regular rewarming and rehydration.
3. Stain with Perls Stain Solution for 15-30min and rinsed with distilled water for 2-5min.(See Note 4)
4. Re-dye with Nuclear Fast Red Solution for 5-10min. Rinse with distilled water for 1-5s.
5. Conventionally dehydrate and transparent, then seal with resinene.

Result

Hemosiderin or Ferric Iron	Blue
Nucleus and Background	Red

Negative Control(Optional)

Take the same control section and dewax it to water, then treat it with a 1% hydrochloric acid solution at room temperature for 20min before staining normally. The result is negative for staining.

Take some ferrous salt and prepare an aqueous solution samples (commonly 0.1% ferrous sulfate aqueous solution), and drop a drop of staining solution, the result is negative for kit.

Take some iron salt and prepare an aqueous solution samples (commonly 0.1% ferric chloride or ferric sulfate aqueous solution). Drop a drop of the staining solution, the result is a positive for kit.

Note

1. Section dewaxing should be as clean as possible. Series of ethanol should be replaced frequently.
2. This staining method for tissue sampling should avoid using acidic fixative or chelating agent treatment that may cause iron ion loss. It is recommended to use G2161-neutral formalin fixative (10%) or P1110-4% tissue cell fixative for tissue fixation.
3. During the whole operation process, the container should be clean and avoid the use of metal iron products. When washing sections and containers, distilled water is suitable, because ordinary water contains iron.
4. This staining requires iron content, and it is recommended to use spleen paraffin sections or iron salts for comparative staining before staining to determine the effectiveness of the reagent.
5. For your health and safety, please wear the experimental clothes and disposable gloves.