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## Prussian Blue Iron Stain Kit (Ferrous Iron, Enhance With DAB) V02

**Cat:** G1429

**Size:** 4×20mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Kit components

Reagent		4×20mL	Storage
Reagent(A):Perls Working Solution	A1:Perls Solution A	10mL	RT, avoid light
	A2:Perls Solution B	10mL	RT
Before use, mix A1 with A2 in equal amount to prepare Perls Working Solution.			
Reagent(B):Incubation Solution		20mL	2-8°C, avoid light
Reagent(C): Enhanced Solution	C1:Enhanced Solution A	1mL	2-8°C, avoid light
	C2:Enhanced Solution B	1mL	2-8°C, avoid light
Before use, mix C1, C2 with 1×PBS in 1:1:18 ratio to prepare Enhanced Working Solution, it is ready to use.			
Reagent(D):Redyeing Solution		20mL	RT, avoid light

### Introduction

The Turnbull's blue staining method is a variant of the Perl's blue staining method, which can display non heme divalent iron in cells and tissues by forming a blue precipitate of Turnbull's blue. However, for non iron overloaded cells and tissues, this method often presents false negatives due to insufficient formation of blue precipitates, which has a good display effect on cell aging or drug-induced cell or tissue iron death.

The Prussian Blue Iron Stain Kit (Ferrous Iron, Enhance With DAB) amplifies and displays the positive signal of ferrous ions using the principle of cascade amplification on the basis of conventional Turnbull's blue staining. It is suitable for ferrous iron staining of samples with less ferrous iron deposition. This product uses high-purity raw materials, with clear and stable dyeing results and high repeatability.

### Self Provided Materials

Thermostat or water bath, wet box, 1 × PBS, distilled water

### Protocols(for reference only)

- Cut the tissue in paraffin sections of 3-7µm. Dewax the slices to water routinely.
- Take out the kit from the refrigerator and rewarm it for 20min to room temperature (25-30 °C), inject water into the wet box and preheat it in a 37 °C incubator for 20min.
- Prepare Perls Working Solution in the ratio of 1:1, drop it onto the slices until the tissues are completely covered, and incubate in a wet box at 37 °C for 20min.
- Take out the slices, gently rinse them with distilled water for three times, 10s each time.
- Drip Incubation Solution onto the slices until the tissue are completely covered, and incubate in a wet box at 37 °C for 10-20min. (see note 4)
- Take out the slices, gently soak with 1×PBS for three times, 60s each time. According to C1: C2:1×PBS=1:1:18 to prepare Enhanced Working Solution, the working solution can be stored at 2-8°C for one week.
- Drip Enhanced Working Solution onto the slices until the tissue are completely covered, and incubate in a wet box at 37 °C for 10-20min. (see note 5)
- Take out the slice, gently soak with 1×PBS for three times, 5s each time, and stain with Redyeing Solution for 3-5min.
- Soak in distilled water for 10min.
- Dehydrate in gradient ethanol, transparent by xylene and seal with resinene.

### Result

Ferrous Ions	Yellow Brown or Brown
Nucleus	Light Blue
Cytoplasm	Light Brown or Colorless

## Note

1. Avoid using acidic fixatives to treat tissues. Chromate treatment will also affect the presence of iron.
2. All working solution are easy to lose effect after a long time. It is recommended to use and prepare them now and use them within 2 hours.
3. The incubation time of Incubation Working Solution should be adjusted appropriately due to tissue size and cell density.
4. To enhance the treatment of working solution, it is recommended to control the staining degree under the microscope, when the positive cells are yellow brown to brown black.
5. If the constant temperature treatment at 37 °C cannot be guaranteed, the treatment time can be appropriately extended or shortened according to the room temperature.
6. For your safety and health, please wear experimental clothes and disposable gloves.

## Related products

- G1420 *Prussian Blue Iron Stain Kit (With Neutral Red)*
- G1422 *Prussian Blue Iron Stain Kit (With Nuclear Fast Red)*
- G1424 *Prussian Blue Iron Stain Kit (With Eosin)*
- G1426 *Prussian Blue Iron Stain Kit (For Cells)*
- G1428 *Prussian Blue Iron Stain Kit (Enhance With DAB)*
- G3310 *Lillie's Ferric Iron Stain Kit*
- G3320 *Lillie's Ferrous Iron Stain Kit*