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Nissl Stain Solution (Toluidine Blue Method)

Cat:G1436

Size:100mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Nissl body is a triangular or elliptical block substance distributed in the cytoplasm of nerve cells. It can be dyed purple-blue by methyl violet, methylene blue, toluidine blue and cresyl violet.

Nissl Stain Solution (Toluidine blue method) has many characteristics: simple operation, stable staining, and wide application. It can be used in paraffin tissue sections of Nissl substance/Nissl bodies, neurons and so on. The existence and disappearance of Nissl is an important indicator of whether nerve cells are damaged. When encephalitis, cerebral ischemia and axonal reaction occur, Nissl will dissolve or even disappear.

Protocol (for reference only)

1. For paraffin section: fix fresh tissue in 10% neutral formalin fixative, then dehydrate and embed, the recommended thickness is 5-10 μ m; For frozen section: fix fresh tissue in 10% neutral formalin fixative, then dehydrate in series of sucrose for 24-72h and embed in OCT embedding reagent, the recommended thickness is 10-15 μ m.
2. For paraffin section: dewax to distilled water. For frozen section: soak the section in distilled water and restore to room temperature.
3. Soak the section into Nissl Stain Solution (Toluidine Blue Method) and incubate at 50-60 $^{\circ}$ C for 20-40min.
4. Slightly wash with distilled water.
5. Differentiate by 95% ethanol rapidly.
6. Absolute ethanol dehydration, transparent by xylene, seal with resinene.

Result

Nissl body	Purple Blue
Nucleus	Light Blue
Background	No Color or Light Blue

Note

1. Nissl dissolves easily in vitro, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain.
2. Tissue fixation plays a very important role. Ethanol or neutral formalin solution can be used for fixation.
3. This staining solution has a good effect on Nissl staining of paraffin tissue sections.
4. Differentiation of 95% ethanol should be carried out rapidly. It is appropriate to observe the slices clearly by naked eye and the background is light blue or colorless.
5. For your safety and health, please wear experimental clothes and disposable gloves.