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Modified Sirius Red Staining Solution (For Collagen Fiber)

Cat: G1473

Size: 100mL/500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Collagen fiber is the most widely distributed and abundant fiber in connective tissue. It is widely distributed in various organs. Type I collagen fibers are mainly bone, skin and tendon; Type II collagen fibers are mainly cartilage collagen; Type III collagen fibers are mainly found in embryonic tissues, adult blood vessels and gastrointestinal tract; Type IV collagen fibers are mainly in the basement membrane. Sirius red and its lining dye solution are strong acid dyes, which are easy to combine with the basic groups in collagen molecules and have strong adsorption. Polarized light microscopy shows that the collagen fibers have the property of positive uniaxial birefringence. When combining with Sirius red, it could enhance the birefringence and improve the resolution, so as to distinguish the two types of collagen fibers.

Modified Sirius Red Staining Solution(For Collagen Fiber) is the Reagent B of G1472 , which can be used in combination with hematoxylin or alone. Mainly used for studying abnormal changes in collagen fibers during various tissue lesions. The ability to dye collagen fibers in tissues such as heart and blood vessels red is helpful for the classification and grading of various fibrotic lesions under polarized light microscopy.

Self Provided Materials

10% formalin fixative solution, ordinary optical microscope or polarized light microscope

Protocol(for reference only)

1. Fix the tissue in 10% formalin fixative solution and embed by routine dehydration.
2. Cut the slices into 3-8 μm . Conventional dewaxing to distilled water.
3. Drop stain with Modified Sirius Red Staining Solution(For Collagen Fiber) for 15-30min. For tissues that are easy to be stained, the dyeing time can be controlled within 5-10min.
4. Rinse slightly with distilled water to remove the staining liquid on the surface of the slice.
5. Dehydrate by series of ethanol and transparent by xylene, finally seal with neutral gum.

Result

1) Optical microscope observations	Collagen fiber	Red
	Muscle fiber	Yellow
2) Polarized light microscope observations	Type I collagen fiber	Strong Orange Yellow or Light Red
	Type III collagen fiber	Green

Note

1. In order to display clearly under the polarized light microscope, the slice thickness should be 6-7 μm .
2. The nucleus without hematoxylin staining may be light red, which is a normal phenomenon, but it will not occur after hematoxylin staining.
3. This method is stable and not easy to fade.
4. For your safety and health, please wear experimental clothes and disposable gloves.