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- [1] Yu Jiang, Dantian Zhu, Wenfeng Liu, et al. Hedgehog pathway inhibition causes primary follicle atresia and decreases female germLine stem cell proliferation capacity or stemness. Stem Cell Research & Therapy. July 2019. (IF 4.627)

Alkaline Phosphatase Stain Kit (Kaplow's/Azo Coupling Method)

Cat: G1480

Size: 4×2mL/4×10mL/4×20mL

Storage: -20°C, avoid light, valid for 6 months.

Kit Components

Reagent		4×2mL	4×10mL	4×20mL	Storage
Reagent(A) : ALP Fixative		2mL	10mL	20mL	RT, avoid light
Reagent(B) :ALP Incubation Solution	B1:AS-BI Solution	1mL	5mL	10mL	-20°C, avoid light
	B2:FBB Solution	1mL	5mL	10mL	-20°C, avoid light
Mix B1:B2 in 1:1 ratio as ALP Incubation Solution.It's ready to use.					
Reagent(C) :Nuclear Fast Red Solution		2mL	10mL	20mL	2-8°C, avoid light
Reagent(D) :Methyl Green Solution		2mL	10mL	20mL	RT, avoid light

Introduction

Alkaline phosphatase (ALP or AKP) is widely distributed in mammalian tissues. The optimum pH for its activity is 9.2-9.8. This enzyme mainly exists in active substance exchange sites (cell membranes):Such as the brush-like margin of intestinal epithelium and proximal convoluted tubule of kidney, the stationary cilia of epididymis epithelium, the capillary cholangium of liver, and the endothelium of arterioles and capillary arteries. It is also found in the endoplasmic reticulum, Golgi complex, swallowing vesicles, lysosomes of intestinal epithelium, neutral granules of neutrophils and smooth muscle cell membranes.

Alkaline Phosphatase Staining Kit (Kaplow Method) is not a metal precipitation method to display alkaline phosphatase activity, but an azo coupling method (also known as simultaneous coupling method).The principle is that under the alkaline condition of pH 9.2-9.8, intracellular alkaline phosphatase can hydrolyze AS-BI phosphate and release phosphoric acid and naphthol, which can form colored products with coupling diazo salts and locate in the cytoplasm.The dye can be used for blood, bone marrow or cell smears, frozen and paraffin sections.Unsettled frozen sections can preserve enzyme activity to the greatest extent possible, while paraffin sections have poor preservation of enzyme activity, and it is not recommended to prepare paraffin sections for staining.The alkaline phosphatase positive site is blue and located in the cytoplasm. This method is more reliable than metal salt precipitation method.

Self Provided Materials

Glass slide, Wet box, Optical microscope

Protocol(for reference only)

Reagent C :Nuclear Fast Red Solution may produce suspended solids or a small amount of precipitation due to flocculation. It is recommended to take supernatant or boil water bath for 5-10min and then air it to 30-40 °C. (see Note 3)

Sample Processing

1. Cell smear: Take fresh whole blood or bone marrow to prepare a smear, place the slide at 30 degrees on the slide, directly in front of the blood or cell droplet, move it back slightly to contact the blood or cell droplet, and make the latter scatter along the lower edge of the slide. Then, slide it smoothly forward along the slide plane until evenly spread, and let it dry naturally.
2. Adherent culture of cells and cell creep: remove the culture medium, wash with 1×PBS 3-4 times, and drain excess water. Remove the cell slides.
3. Frozen slice: Remove from the refrigerator and place at room temperature. Soak in distilled water and wash for 1-2 min, then drain excess water.

Staining Incubation

1. Fix with pre-cooled ALP Fixative at 2-8 °C for 30s-3min, and in most cases, 30-60s is sufficient.
2. Wash with distilled water and drain excess water (not too dry).
3. (Slicing treatment) Drip ALP incubation solution onto the slices and incubate them in a wet box at 37 °C for 15-20 min. Wash with distilled water.
4. (Adherent cells) Drip ALP incubation solution to cover the cells, incubate at 37 °C for 15-20 min, and wash with distilled water.
5. Re staining: Stain with Nuclear Fast Red Solution or Methyl Green Solution for 3-5 min. For frozen section,

1-2min is sufficient.

Sealing observation

1. (Slicing treatment) Wash slightly with distilled water, drain excess water, add water-based sealing agent dropwise, seal the film, and observe under a microscope.
2. (Adherent cells) Wash slightly with distilled water and observe under a microscope with water or 1×PBS.

Result

Positive site of ALP	Blue
Nucleus	Red or Green

Criteria for blood and bone marrow smears judging

Judge according to 100 neutrophil-positive granules and report the level in 0-4

Positive level	Dyeing Characteristics
0	No Granule
1	Slightly Granule
2	Medium-Grade Granule
3	Majority Granule
4	Full Granule

Note

1. Blood or bone marrow smears or other samples should be fresh. Cut sections in appropriate thickness and fix it in time to avoid reducing activity of the enzyme.
2. For cultured cells and cell slides, if there is cell floating phenomenon after fixation, 4% paraformaldehyde can be used instead of ALP fixation solution for fixation treatment, and fixed at room temperature for 10-15 minutes.
3. Reagent C :Nuclear Fast Red Solution is a colloidal solution, which is stored at low temperature (lower than 25 °C) or stored for a long time. Suspended solids or a small amount of precipitation are generated due to flocculation, which is a normal phenomenon and generally does not affect the use. If the colloid solution is evenly dispersed in the boiling bath, tighten the bottle cap for 5-10min to recover the turbid solution.
4. In the process of cultured cells staining, slightly wash and stain to avoid damage or loss of cells.
5. ALP Incubation Solution is easy to lose effect. Use it as soon as possible.
6. Re-dyeing by Nuclear Fast Red Solution or Methyl Green Solution.
7. There should be a positive control for each sample.

Reference

[1] Yu Jiang, Dantian Zhu, Wenfeng Liu, et al. Hedgehog pathway inhibition causes primary follicle atresia and decreases female germLine stem cell proliferation capacity or stemness. Stem Cell Research & Therapy. July 2019. (IF 4.627)