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Cat: G1491

Size: 4×2mL/4×10mL/4×20mL

Storage: -20°C, avoid light, valid for 6 months.

Kit Components

Reagent		4×2mL	4×10mL	4×20mL	Storage
Reagent(A): ACP Fixative		10mL	50mL	100mL	2-8°C, avoid light
Reagent(B): ACP Incubation Solution	B1:AS-BI Buffer	0.2mL	1mL	2×1mL	-20°C, avoid light
	B2: GBC Solution	0.02mL	0.1mL	0.2mL	-20°C, avoid light
	B3: ACP Buffer	1.8mL	9mL	18mL	RT, avoid light
Mix reagent B1, B2, B3 in 10:1:90 ratio as ACP Incubation Solution before use.					
Reagent(C): Hematoxylin Solution		2mL	10mL	20mL	2-8°C, avoid light
Reagent(D): Methyl Green Solution		2mL	10mL	20mL	RT, avoid light
Reagents C and D are both re staining solutions, and one can be used during the staining process. Repeated staining is not recommended.					

Introduction

Acid phosphatase (ACP) is widely distributed in various tissues, mainly in the lysosome, so it is often used as a lysosome marker enzyme. The acid phosphatase outside of lysosome exists in endoplasmic reticulum and cytoplasm. Acidic phosphatase is different in all kinds of animals. The optimum pH for its activity is 4.5-5.5.

The reaction principle of Acid Phosphatase Stain Kit (Kaplow's/Azo Coupling Method) is that phosphoric acid and naphthol are released by hydrolysis of acid phosphatase with AS-BI as substrate at acidic pH. Naphthol was coupled with diazo salts to form colored products, which were localized in the cytoplasm. It usually used for fresh blood smear, cell smear, frozen section, etc.

Protocol (for reference only)

Sample Processing

1. Cell smear: Take fresh whole blood or bone marrow to prepare a smear, place the slide at 30 degrees on the slide, directly in front of the blood or cell droplet, move it back slightly to contact the blood or cell droplet, and make the latter scatter along the lower edge of the slide. Then, slide it smoothly forward along the slide plane until evenly spread, and let it dry naturally.
2. Adherent culture of cells and cell creep: remove the culture medium, wash with 1×PBS 3-4 times, and drain excess water. Remove the cell slides.
3. Frozen slice: Remove from the refrigerator and place at room temperature. Soak in distilled water and wash for 1-2 min, then drain excess water.
4. Paraffin section: Paraffin sectioning is dewaxed for 5-10 min and repeated once. Anhydrous ethanol for 5 min, 90% ethanol and 70% ethanol for 2 min each, soak in distilled water for 2 min, and drain excess water.

Staining Incubation

1. Fix with pre-cooled ACP Fixative at 2-8 °C for 30s-3min, and in most cases, 30-60s is sufficient.
2. Wash with distilled water and drain excess water (not too dry).
3. (Slicing treatment) Drip ACP incubation solution onto the slices and incubate them in a wet box at 37 °C for 45-60 min. Wash with distilled water.
4. (Adherent cells) Drip ACP incubation solution to cover the cells, incubate at 37 °C for 45-60 min, and wash with distilled water.
5. Re staining: Stain with Hematoxylin Solution or Methyl Green Solution for 2-3 min. If Hematoxylin Solution is used for re staining, must return blue with tap water for 10 min or with 1×PBS for 3-5 min after staining.

Sealing observation

1. (Slicing treatment) Wash slightly with distilled water, drain excess water, add water-based sealing agent dropwise, seal the film, and observe under a microscope.
2. (Adherent cells) Wash slightly with distilled water and observe under a microscope with water or 1×PBS.

Result

Positive Site	Purplish Red
Nucleus	Blue(Hematoxylin) or Green(Methyl Green)

Clinical Significance

1. ACP staining of hairy cells in hairy cell leukemia is strongly or moderately positive, and is not inhibited by tartaric acid.
2. ACP staining is positive in immature monocytes of acute leukemia, and weak positive in polymorphocytes. The response of progranulocytes to ACP is different.
3. Tlymphocyte ACP staining is positive, with large and densely distributed granules. B lymphocyte is negative or weak positive with small granules.
4. Gaucher cells are strongly positive and Niemann-Pick cells are negative or weak positive.

Note

1. This staining kit is an enzyme active in situ staining kit. The reaction substrate and dye are prone to failure, and it is recommended to store the product properly according to the temperature described in the instruction manual as soon as received.
2. When staining frozen sections, the exposure time of sections at room temperature should be reduced.
3. Enzyme activity is easily affected by various factors, leading to attenuation or disappearance. It is recommended to take fresh samples and use pre cooled fixing solution to fix them in a refrigerator at 2-8 °C for no more than 24 hours. It is recommended to use low melting point wax (52-56 °C) during embedding to avoid thermal deactivation.
4. For your safety and health, please wear laboratory clothes and disposable gloves for operation.