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Cat: G1511

Size: 3×50mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×50mL	Storage
Reagent(A): Mayer Hematoxylin Staining Solution	50mL	RT, avoid light
Reagent(B): Sudan III Differentiation Solution	50mL	RT
Reagent(C): Sudan III Solution	50mL	RT, avoid light

Introduction

Neutral fat is a lipid composed of three molecules of fatty acids and one molecule of glycerol, which is neutral. Neutral fats are one of the ways to store energy, releasing energy during oxidation. Under normal circumstances, except for adipocytes, other cells can hardly see lipid droplets under an optical microscope. If a large number of lipid droplets appear in the cytoplasm, it is considered steatosis, commonly seen in liver cells, myocardial cells, renal tubular epithelial cells, etc. Neutral fat staining often uses methods such as Sudan II, Sudan III, Sudan IV, Sudan Black B, and Oil Red O. The mechanism of Sudan dye lipid staining is generally believed to be purely physical, involving lipid solubility and adsorption. Sudan dyes, due to their higher solubility in lipids than in organic solvents, transfer the dye from the dye solution to the dyed lipids during dyeing, resulting in the lipids presenting the color of the dye solution.

Sudan III Stain Kit is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The sample should not fixed with ethanol (10% formalin if required) or paraffin. Frozen sections or carbon wax sections were needed.

Self Provided Materials

Slide, Microscope, 70% alcohol, Distilled Water, S2150-Glycerol Gelatin or S2155-Aqueous Mounting Medium(PVP)

Protocol (for reference only)

1. Prepare low-temperature sections of fresh tissue, generally at -20°C to -25°C. If the sample is lipoma, it should be adjusted to -30 °C.
2. Cut the frozen section at 6-15um (6-8 um is preferable) and attach to the slide. Wash with 70% alcohol.
3. Stain with Sudan III Solution for 10min. Wash with 70% alcohol to remove excess solution. Wash with distilled water for 1 min.
4. Re-dyeing slightly with Mayer Hematoxylin Staining Solution for 1-2min.
5. (optional) Under the view of microscope, if the staining is too deep, can differentiate by Sudan III Differentiation Solution for several seconds. Rinse with tap water for 5-10 min to bluing.
6. Use filter paper to absorb the slices and surrounding water and let them dry slightly.
7. Glycerol gelatin or S2155-Aqueous Mounting Medium(PVP) sealing.

Result

Neutral fat	Orange Red
Nucleus	Light Blue

Note

1. Samples should not be fixed with ethanol or embedded with paraffin, and frozen sections should be used. It is recommended to use 10% neutral formalin or 10% formaldehyde-calcium solution to fix.
2. Precipitation of dyes must be prevented during dyeing. Therefore, slices should be covered when they are put into dye solution, so as not to contact with flowing air to avoid precipitation when the solution volatilizes.
3. Frozen sections are easier to be colored, and over-dyeing should be avoided when re-dyeing with Mayer Hematoxylin Staining Solution.
4. Samples sealed with glycerol gelatin did not last long. If long-term preservation is required, the edge of the boundary between the cover slide and the slide can be sealed with neutral gum.
5. For your safety and health, please wear lab clothes and disposable gloves.