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**Cat:** G1540

**Size:** 3×100mL/3×500mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	3×100mL	3×500mL	Storage
Reagent(A): Diff-Quik Fixative	100mL	500mL	RT, avoid light
Reagent(B): Diff-QuikI	100mL	500mL	RT, avoid light
Reagent(C): Diff-QuikII	100mL	500mL	RT, avoid light

### Introduction

Diff-Quik Staining is a fast staining method based on Wright Staining, which is one of the common staining methods in cytology. The dye solution is prepared by the rapid staining method recommended by the World Health Organization (WHO). Similar to Wright stain, it is improved by Romanowsky Stain technology principle. The staining result is also very similar to that of Wright's, but the time required for Diff Quik Staining is very short, and the staining can be completed within 90s.

Diff-Quick Stain is mainly used for blood cell smear, bone marrow smear, vaginal secretion smear and exfoliated cell smear. Diff-Quick Stain Kit (With Fixative) is very suitable for batch dyeing with clear background and no sediment.

### Protocol(for reference only)

1. Prepare blood smear or bone marrow smear by routine method, and fix in Diff-Quik Fixative for 20s after natural drying.
2. Staining with Diff-Quik I for 5-10s (lifting the section up and down for 2-3 times to evenly distribute the dye solution), take it out immediately without washing.
3. Staining with Diff-Quik II for 10-20s (lifting the section up and down for 2-3 times to evenly distribute the dye solution), take it out immediately.
4. After washing with distilled water, quickly observe under the microscope immediately when it is wet, or observe after dehydration in 95% ethanol for 5s, 100% ethanol (I) for 5s, 100% ethanol (II) for 30s, transparency by xylene for 1min and sealing.

### Result

Nucleus, Leukocyte	Deep Blue
Matrix, Lymphocyte	Purple
Cytoplasm, Red Blood Cell	Pink Red

### Note

1. Diff-Quik Fixative is a specific fixative for Diff-Quik. If it is used in large quantity, can replace by methanol.
2. Blood smear or bone marrow smear should be uniform in thickness, so as not to affect the staining effect. Blood cell smear staining requires fresh whole blood or EDTA anticoagulant.
3. Bone marrow smear staining requires the prepared smear to be shaken or dried quickly in the air, so as to prevent cell shrinkage and deformation or hemolysis due to air humidity. It is not allowed to dry the smear with high temperature or fire baking.
4. Smear staining of vaginal secretions requires fresh specimens to be fixed with flame or alcohol as soon as possible after smearing to avoid cell deformation.
5. Please do not remove the dye solution or wash the smear directly. Do not pour out the dye first, so as to prevent the dye from settling on the smear.
6. The dye solution can be used for limited reuse. If there is sediment, should filter it before use.
7. If the dye is too deep, it can be decolorized properly with methanol or alcohol, and it is better not to be re dyed.
8. If the dyeing is too deep or too shallow, adjust the dyeing time or working solution concentration.
9. The pH value has certain influence on dyeing. The section should be clean and free of acid and alkali pollution, so as not to affect the dyeing effect.