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
Collagen Fiber And Elastic Fiber Staining Kit(EVG-Modified Weigert Method)^{V02}

Cat: G1591

Size: 5×50mL/5×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent		5×50mL	5×100mL	Storage
Reagent(A): Weigert Oxidant	A1: Weigert Oxidant A	25mL	50mL	RT, avoid light
	A2: Weigert Oxidant B	25mL	50mL	RT
Mix A1 with A2 in equal amount to form Weigert Oxidant. It is recommended to use it within 3 hours.				
Reagent(B): Weigert Bleaching Solution		50mL	100mL	RT
Reagent(C): Weigert Fuchsin Staining Solution		50mL	100mL	2-8°C, avoid light
Reagent(D): Acid Differentiation		50mL	100mL	RT
Reagent(E): Modified VG Staining Solution		50mL	100mL	RT, avoid light

Introduction

Elastic fibers are mainly distributed in the arterial walls, alveolar walls, and skin of the human body. When fresh, they are yellow and have strong refractive properties. In routine tissue staining, it is often necessary to simultaneously display two types of fibers to study the coexistence relationship between fibers and their dependence on other structures. Based on this, an EVG composite staining kit was developed on the basis of a separate fiber staining kit, which is used in combination with conventional elastic fiber staining and VG collagen staining.

This reagent kit adopts a combined staining mechanism of weigert fuchsin and modified VG staining solution, which can simultaneously stain elastic fibers and collagen fibers on a single slice without fading easily. After staining, the elastic fibers are purple black, the collagen fibers are red, and the muscles are yellow. It is suitable for observing whether there are abnormal lesions in the elastic fibers and collagen fibers in the tissue, thus assisting in diagnosis. This reagent kit has excellent quality, simple operation, stable staining, and beautiful and clear colors.

Protocol(for reference only)

1. Take fresh samples, fix and embed with paraffin. Cut the section into 3-8um.
2. After 50min of heat treatment on paraffin slices, dewax by xylene for 20min, then treat with gradient ethanol to distilled water for each step 5min.
3. Treat with Weigert Oxidant for 3min, then wash with distilled water for 1min.(see note 2)
4. Treat with Weigert Bleaching Solution for 3min, then fully wash with distilled water for 1min.(see note 2)
5. Soak the section into Weigert Fuchsin Staining Solution and cover at room temperature for 1-3 hours, or stain at 56 °C for 1 hour.
6. Differentiate with Acid Differentiation until no dye is removed. Then wash with distilled water.
7. Dye with Modified VG Staining Solution for 5-10min, and wash with distilled water for 30s.
8. 75%, 85%, 95%, and 100% ethanol were subjected to conventional dehydration for 3-5 seconds each. Xylene was transparent twice for 1min each time, and neutral gum was used for sealing.
9. Microscopic observation of staining results.

Result

Elastic fiber	Brown or Purple Black
Collagen fiber	Red
Muscle fiber, Red blood and Background	Yellow

Note

1. Most fixing solutions are suitable, but fixing solutions containing chromium salts can result in lighter staining and easier dispersion. Slice dewaxing should be as clean as possible, and series of ethanol should be frequently replaced with a new solution.
2. If you only want to display mature elastic fibers, you can skip the oxidation bleaching step and proceed with

subsequent dyeing directly.

3. If preheating (56 °C) is not added before dyeing Weigert Fuchsin Staining Solution, the dyeing time can be appropriately extended
4. The Acid Differentiation usually takes only 3-10 seconds to remove the background color, and the differentiation must be terminated as soon as possible when the color is clear through visual inspection.
5. The Modified VG Staining Solution is stable in dyeing, not easy to fade after washing, and can be thoroughly cleaned.
6. For your safety and health, please wear laboratory clothes and disposable gloves for operation.