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## Elastic Fiber Stain Kit(Aldehyde-Fuchsin Method)

**Cat:** G1593

**Size:** 4×50mL

**Storage:** 2-8°C, avoid light, valid for 3 months.

### Kit Components

Reagent		4×50mL	Storage
Reagent(A): Acid Oxidizing Solution	A1: Oxidizing Solution	25mL	2-8°C, avoid light
	A2: Acidizing fluid	25mL	2-8°C
Mix reagent A1, A2 in ratio 1:1 to form Acid Oxidizing Solution before use. Please use it in time.			
Reagent(B): Acid Bleach Solution		50mL	RT
Reagent(C): Aldehyde-Fuchsin Solution		50mL	2-8°C, avoid light
Reagent(D): Orange G Solution		50mL	RT, avoid light

### Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elasticize usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin- fuchsin , and Weigert's elastic stain in histological sections.

The dyeing principle of Elastic Fiber Stain Kit(Aldehyde-Fuchsin Method) is that mature fuchsin aldehyde has a strong affinity for special proteins such as elastic fibers and mucopolysaccharides containing sulfate. The dye solution can be used to display elastic fibers, mast cell particles, lipofuscin, eosinophils, etc.

### Protocol(for reference only)

1. Fix in 10% neutral formalin fixative, and routinely dehydrate and embed.
2. Cut the paraffin embedded tissue into 4µm thin sections and routine dewax to water.
3. Rinse the section with Acid Oxidizing Solution just prepared for 5min. Wash with distilled water for 1-2min.
4. Bleaching with Acid Bleach Solution for 5min. Wash with distilled water for 1-2 min.
5. Rinse with 70% alcohol.
6. Soak the section in Aldehyde-Fuchsin Solution and cover dyeing at RT for 1-2h or at 56°C for 15-30min.
7. Quickly rinse with 70% alcohol for 3-5s and wash with distilled water to stop differentiation until no more solution.
8. Stain with Orange G Solution for 2-3s. Wash with distilled water for 1-2min.
9. Quickly dehydrate with 95% alcohol and absolute alcohol, transparent with xylene and seal.

### Result

Elastic Fiber	Purple to dark purple
Mast cell granules and mucus	Purple to dark purple
background	Different degrees of yellow

### Note

1. This method can display elastic fiber, anterior elastic fiber and acid resistant fiber, but it needs to slice slightly thicker, 7µm is the best.
2. The oxidizing solution and acidifying solution should not be mixed in advance, and it is better to use them immediately.
3. When aldehyde fuchsin is dyed, it should be capped to prevent the solution from volatilizing.
4. When the Aldehyde-Fuchsin Stain Solution is stored for a long time, the dyeing power will decrease, and the dyeing time should be increased.
5. When islet beta-cells were stained, the time should be controlled at 30 min; when pituitary basophils were stained, the time is controlled at 60 min.
6. Orange G dyeing should be light, otherwise it will cover up the color of elastic fiber.
7. For your safety and health, please wear experimental clothes and disposable gloves.