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Collagen Fiber And Elastic Fiber Staining Kit(EVG-Victoria Method)

Cat: G1595

Size: 4×50mL/4×100mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent		4×50mL	4×100mL	Storage
Reagent(A):Tanake Oxidant	A1: Tanake Oxidant A	25mL	50mL	RT, avoid light
	A2: Tanake Oxidant B	25mL	50mL	RT
Before use, mix A1 with A2 in equal amount to form Tanake Oxidant, it is not suggested to prepare in advance.				
Reagent(B): Tanake Bleaching Solution		50mL	100mL	RT
Reagent(C): Victoria Blue Staining Solution		50mL	100mL	2-8°C, avoid light
Reagent(D): Modified VG Staining Solution		50mL	100mL	RT, avoid light

Introduction

Elastic fibers are mainly distributed in the arterial walls, alveolar walls, and skin of the human body. When fresh, they are yellow and have strong refractive properties. In routine tissue staining, it is often necessary to simultaneously display two types of fibers to study the coexistence relationship between fibers and their dependence on other structures. Based on this, an EVG composite stain was developed on the basis of a separate fiber staining kit, which is used in combination with conventional elastic fiber staining and VG collagen staining.

This reagent kit adopts a combined staining mechanism of victoria blue staining solution and modified VG staining solution, which can simultaneously stain elastic fibers and collagen fibers on a single slice without fading easily. After staining, the elastic fibers are blue, the collagen fibers are red, and the muscles are yellow. It is suitable for observing whether there are abnormal lesions in the elastic fibers and collagen fibers in the tissue, thus assisting in diagnosis. This reagent kit has excellent quality, simple operation, stable staining, and beautiful and clear colors.

Protocol(for reference only)

1. Take fresh samples, fix and embed with paraffin. Cut the section into 3-8um.
2. After 50min of heat treatment on paraffin slices, dewax by xylene for 20min, then treat with gradient ethanol to distilled water for each step 5min.
3. Treat with Tanake Oxidant for 3min, then wash with distilled water for 1min.
4. Treat with Tanake Bleaching Solution for 3min, then fully wash with distilled water for 1min.
5. Slightly wash with 70% ethanol. Soak the section into Victoria Blue Staining Solution and cover dyeing at room temperature for 15-30min.(see note 3)
6. Soak the slices in 70% ethanol twice, each time for 10 seconds. During this time, lift and place the slices down until no staining solution is removed from the slices, and wash them slightly with running water..
7. Dye with Modified VG Staining Solution for 5-10min, and wash with distilled water for 30s.
8. 75%, 85%, 95%, and 100% ethanol were subjected to conventional dehydration for 3-5 seconds each. Xylene was transparent twice for 1min each time, and neutral gum was used for sealing.
9. Microscopic observation of staining results.

Result

Elastic fiber	Blue
Collagen fiber	Red
Muscle fiber, Red blood and Background	Yellow or Yellowish Green

Note

1. Slice dewaxing should be clean, and series ethanol should be frequently replaced with a new solution.
2. If you only want to display mature elastic fibers, you can skip the oxidation bleaching step and proceed with subsequent dyeing directly.
3. Victoria Blue Staining Solution is stable, stored for a long time, and can be reused.
4. The active ingredients in Victoria Blue Staining Solution are easily volatile, and immersion dyeing is preferred. The immersion time can be appropriately extended according to the actual dyeing depth required.
5. The Modified VG Staining Solution is stable in dyeing, can be thoroughly cleaned.
6. For your safety and health, please wear laboratory clothes and disposable gloves for operation.