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## Elastic Fiber Stain Kit(Victoria Blue Method)

**Cat:** G1596

**Size:** 4×50mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent		4×50mL	Storage
Reagent(A): Tanake Oxidizing Solution	A1: Tanake Oxidizing Solution A	25mL	2-8°C, avoid light
	A2: Tanake Oxidizing Solution B	25mL	RT
Before use, mix reagent A2, A1 in ratio 1:1 to form tanake oxidant, which should be used in time.			
Reagent(B): Tanake Bleach Solution		50mL	RT
Reagent(C): Victoria Blue Solution		50mL	RT, avoid light
Reagent(D): Nuclear Fast Red Solution		50mL	RT, avoid light

### Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elasticized usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin- fuchsin , and Weigert's elastic stain in histological sections. Victoria blue is a kind of phenylmethane dye, which can show elastane and hepatitis B surface antigen. It is not easy to over dye and fade after dyeing.

Elastic Fiber Stain Kit(Victoria Blue Method) is based on the improvement of tanake Victoria blue method. It is applicable to display elastic fiber and hepatitis B surface antigen. It is used to observe whether the elastic fiber in the tissue proliferates or breaks and disintegrates clinically, so as to assist in diagnosis.

### Protocol (for reference only)

*Reagent(D): Nuclear Fast Red Solution may produce suspended solids or a small amount of precipitation due to flocculation. It is recommended to take supernatant or boil water bath for 5-10min and then air it to 30-40 °C. (see Note 2)*

1. Fix in 10% neutral formalin fixative, and routinely dehydrate and embed.
2. Cut the paraffin embedded tissue into 4μm thin sections and routine dewax to water.
3. Treat the section with Tanake Oxidizing Solution for 5min. Wash with distilled water for 1min.
4. Treat the section with Tanake Bleach Solution for 2min. Wash with distilled water for 1 min.
5. Rinse with 70% alcohol. Stain the section in Victoria Blue Solution with cap for 24h.
6. Rinse with 70% alcohol twice till no stain solution comes out. Wash with distilled water for 1min.
7. Stain with Nuclear Fast Red Solution for 5-10min. Wash with distilled water for 1min.
8. Dehydrate with series alcohol, transparent with xylene and seal.

### Result

Elastic Fiber	Blue
Nucleus	Red

### Note

1. Most fixatives can be used. However, when the fixative used containing chromium salt, the staining is light and easy to disperse.
2. Reagent(D): Nuclear Fast Red Solution is a colloidal solution, which is stored at low temperature (lower than 25 °C) or stored for a long time. Suspended solids or a small amount of precipitation are generated due to flocculation, which is a normal phenomenon and generally does not affect the use. If the colloid solution is evenly dispersed in the boiling bath, tighten the bottle cap for 5-10min to recover the turbid solution.
3. Victoria blue staining solution is stable and can be used repeatedly.
4. The staining time of Victoria blue could be from 24 h to 48 h, most of which is 24 h.
5. For your safety and health, please wear experimental clothes and disposable gloves.