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规格: 100mL/500mL

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Cat: G1611

Size: 100mL/500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Papanicolaou method is widely used in routine cytological staining. Initially, Papanicolaou Stain was only used to detect estrogen levels in the vaginal epithelium and pathogens such as Candida and Trichomonas in the reproductive tract. Its principle was to use orange G6 in combination with EA36 or EA50 to dye the cytoplasm into bright green, blue, and pink colors. At present, the improved Pap staining solution contains multiple ions and has polychromatic staining efficiency. After staining, the cytoplasm is bright, transparent, and the nuclear membrane, nucleolus, and chromatin structure are clear. The nuclear staining solution is mainly Harris hematoxylin staining solution, while the cytoplasmic staining solution is mainly EA36 staining solution, EA50 staining solution, and EA65 staining solution. Pap staining solution can also be used for cell exfoliation specimens. After staining, the nucleus appears blue or black, the cytoplasm of keratinized squamous cells appears pink or orange red, and the cytoplasm of non keratinized cells appears green or blue-green.

Papanicolaou EA50 Stain Solution can make the cytoplasm of different degrees of keratinization have different colors. It is not only suitable for gynecological cytological smear staining, such as screening cervical cancer and precancerous lesions, but also suitable for the staining of pleural fluid, ascites, sputum and other non gynecological cell samples. It is usually used with hematoxylin.

Self Provided Materials

Fixative(like AF Fixative) ,Series of alcohol, Microscope, Acid alcohol differentiation solution

Protocol(for reference only)

1. For cell smear, fix with 95% alcohol for 10-15mins.
2. (optional)Rinse in 75% alcohol for 1min. Rinse with tap or distilled water for 1min.
3. (optional)Stain with Hematoxylin Solution for 5-10mins.Rinse with tap water for 2mins.
4. (optional)Differentiate with 1% acid alcohol differentiation solution for about 4-5s or 0.5% acid alcohol differentiation solution for about 10s.Rinse with tap water for 2mins.
5. (optional)Blue with Bluing Solution for 2mins.Rinse with tap water for 2mins.
6. (optional)Add 95% ethanol dropwise to cover the smear and balance for 5-10s.
7. (optional)Stain with Orange G Solution for 2mins. Rinse with 95% ethanol for 5-10s to remove excess dye solution.
8. Stain with EA50 Solution for 3-5mins. Rinse with 95% ethanol for 5-10s to remove excess dye solution.
9. Dehydrate in absolute alcohol for 1min twice.
10. Transparent with xylene and seal with resinene.

Result

Nucleus	Blue Purple or Black
Cytoplasm of non keratinized cells	Light Blue or Light Green
Cytoplasm of keratinocytes	Pink or Orange

Note

1. The dye solution shall be filtered before use. If it is reused, the dye solution shall be changed frequently.
2. This product contains a certain proportion of alcohol components. Pay attention to cover and stain to prevent volatilization, and tighten the bottle neck in a timely manner when taking it.
3. For your safety and health, please wear experimental clothes and disposable gloves.