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规格: 4×100mL/4×500mL

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Cat: G1612

Size: 4×100mL/4×500mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	4×100mL	4×500mL	Storage
Reagent(A): Hematoxylin Solution	100mL	500mL	RT, avoid light
Reagent(B): Bluing Solution	100mL	500mL	RT
Reagent(C): Orange G Solution	100mL	500mL	RT, avoid light
Reagent(D): EA36 Solution	100mL	500mL	RT, avoid light

Introduction

Papanicolaou method is widely used in routine cytological staining. Initially, Papanicolaou Stain was only used to detect estrogen levels in the vaginal epithelium and pathogens such as Candida and Trichomonas in the reproductive tract. Its principle was to use orange G6 in combination with EA36 or EA50 to dye the cytoplasm into bright green, blue, and pink colors. At present, the improved Pap staining solution contains multiple ions and has polychromatic staining efficiency. After staining, the cytoplasm is bright, transparent, and the nuclear membrane, nucleolus, and chromatin structure are clear. The nuclear staining solution is mainly Harris hematoxylin staining solution, while the cytoplasmic staining solution is mainly EA36 staining solution, EA50 staining solution, and EA65 staining solution. Pap staining solution can also be used for cell exfoliation specimens. After staining, the nucleus appears blue or black, the cytoplasm of keratinized squamous cells appears pink or orange red, and the cytoplasm of non keratinized cells appears green or blue-green.

Papanicolaou EA36 Stain Kit use EA36 staining solution for cytoplasmic staining and non-toxic modified hematoxylin staining solution for nuclear staining. It is more suitable than EA50 for gynecological cytology smear staining, such as screening cervical cancer and precancerous lesions.

Self Provided Materials

Fixative(like AF Fixative) ,Series of alcohol, Microscope, Acid alcohol differentiation solution

Protocol(for reference only)

1. For cell smear, fix with 95% alcohol for 10-15mins.
2. Rinse in 75% alcohol for 1min. Rinse with distilled water for 1min.
3. Stain with Hematoxylin Solution for 5-10mins. Rinse with distilled water for 1min.
4. (Optional) Differentiate with 1% acid alcohol differentiation solution for about 4-5s or 0.5% acid alcohol differentiation solution for about 10s. Rinse with tap water for 1min.
5. Blue with Bluing Solution for 3-5mins. Rinse with tap water for 1min.
6. Add 95% ethanol dropwise to cover the smear and balance for 5-10 seconds.
7. Stain with Orange G Solution for 2mins. Rinse with 95% ethanol for 5-10 seconds to remove excess dye solution. (See Note 3)
8. Stain with EA36 Solution for 3-5mins. Rinse with 95% ethanol for 5-10 seconds to remove excess dye solution. (See Note 3)
9. Dehydrate in absolute alcohol for 1min twice.
10. Transparent with xylene and seal with resinene.

Result

Nucleus	Blue Purple or Black
Cytoplasm of non keratinized cells	Light Blue or Light Green
Cytoplasm of keratinocytes	Pink or Orange

Note

1. All dye solutions should be filtered before use and changed frequently.
2. If nuclear staining is not performed, steps 3, 4, and 5 can be omitted, and the rehydration process cannot be omitted.
3. Steps 7 and 8 should avoid dry tablets, otherwise sedimentation may occur and affect observation.
4. For your safety and health, please wear experimental clothes and disposable gloves.