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Cat: G2030

Size: 100mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Melanin is a kind of Nonhematogenous endogenous pigment, which is a group of pigments from light brown to black. This pigment is usually found in the skin, eyes, the substantia nigra of the brain and hair follicles. Lipofuscin is a granular brown yellow pigment, which is composed of the remains containing fat and lysosomal digests. It is believed that lipofuscin is produced by the oxidation of lipids and lipoproteins. The oxidation process of lipofuscin is slow and gradual, so the pigment will show different dyeing reactions, different colors, and different shapes and sizes.

Melanin And Lipofuscin Stain Solution(Nile Blue Method) can be used to display melanin and lipofuscin by Nile blue. This method is easy to operate, but this solution may not work well for some samples.

Protocol(for reference only)

1. Fix with 10% neutral formalin, routinely dehydrate and embed, cut into slices. The thickness of the slice is recommended to be 3-5um.
2. Deparaffinize experimental and control sections to distilled water.
3. Stain with Melanin and Lipofuscin Stain Solution (Nile Blue Method) for 20-30mins.
4. Wash with distilled water three times and each time for 1-2min.
5. Seal with glycerin gelatin.

Result

Melanin	Dark Green
Lipofuscin	Dark Blue or Turquoise
Nucleus	Blue or Colorless

Note

1. This solution may not work well for some samples.
2. When staining frozen sections, neutral lipids such as triglycerides, cholesterol esters and steroids are dyed red or pink, and acid fats such as fatty acids and phospholipids are dyed blue.
3. For your safety and health, please wear experimental clothes and disposable gloves.