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Cat: G2032

Size: 3×50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	3×50mL	Storage
Reagent(A): Fontana Ammoniac Silver Solution	50mL	2-8°C, avoid light
Reagent(B): Hypo Solution	50mL	RT
Reagent(C): Neutral Red Staining Solution	50mL	RT, avoid light

Introduction

Melanin is a kind of Nonhematogenous endogenous pigment, which is a group of pigments from light brown to black. This pigment is usually found in the skin, eyes, the substantia nigra of the brain and hair follicles. Melanin has a significant physical property. Because the melanin formed in melanosome is closely bound with protein, melanin is completely insoluble in most organic solvents. Another physical property of melanin is that it can be slowly bleached by strong oxidant. Under pathological condition, this kind of pigment can also appear in benign nevus cell tumor and malignant melanoma. Many methods can be used to identify melanin and melanogenesis cells, such as reduction method: Masson Fontana silver method and Schmorl trivalent iron potassium ferricyanide reduction experiment, enzyme method (such as DOPA reaction), fluorescence method, immunohistochemistry, etc.

Masson-Fontana Melanin Stain Kit shows melanin by the argentaffin reaction, that is, melanin has the property of reducing silver ammonia solution to metallic silver. After dyeing, melanin is black. This method belongs to the commonly used special melanin dyeing method, and the effect is better than other dyeing methods.

Self Provided Materials

10% Neutral Formalin Fixative, Incubator or Water Bath, Distilled Water

Protocol(for reference only)

1. Fix with 10% neutral formalin, routinely dehydrate and embed, cut into slices. The thickness of the slice is recommended to be 3-5um.
2. Deparaffinize experimental and control sections to distilled water.
3. Soak in Fontana Ammoniac Silver Solution avoid light at room temperature for 12-24h or 56°C incubator for 30-40min.
4. Wash with distilled water for five-six times fifth and each time for 1-2min.
5. Treat with Hypo Solution for 1-5min. Wash with distilled water for 3-5min.
6. (Optional)Re-dye with Neutral Red Staining Solution for 5mins. Wash with distilled water for 1min.
7. Conventional dehydration and transparency, seal with resinene.

Result

Melanin, Argyrophil Granules	Black
Nucleus	Red

Note

1. Pay attention to control the staining time of Fontana Ammonia Silver Solution under the microscope, and it is better for melanin particles to appear black. If the staining time is too long, lipofuscin, bile pigment and orange blood will also be black, which is not easy to distinguish melanin.
2. Utensils in contact with Fontana Ammoniacal Silver Solution should be cleaned strictly, otherwise there may be no chemical reaction of residual pollutants.
3. Samples exposed to 56 °C for a long time may produce fine precipitates.
4. The fragile sample may be coated with a layer of collodion to prevent the slice from lifting from the slide.
5. For your safety and health, please wear experimental clothes and disposable gloves.