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**Cat:** G2041

**Size:** 3×50mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit components

Reagent	3×50mL	Storage
Reagent(A): Gram Eosin Staining Solution	50mL	RT, avoid light
Reagent(B): Methyl Violet Staining Solution	50mL	RT, avoid light
Reagent(C): Gram Iodine Solution	50mL	RT, avoid light

### Introduction

Pathological endogenous precipitates are a part of pigmentation. After certain pathological changes, tissue cells form precipitates with different shapes and characteristics. The special protein formed by this polymerization can show cellulose protein after staining. Cellulose is a special protein formed by the polymerization of fibrin molecules in the blood, also known as fibrin. This protein exists in the tissue in the form of curved filament cellulose, most of which are in a network structure, sometimes in a thick cellulose network, old agglutinable and amorphous blocks.

Fibrin Stain Kit (Gram Method) is mainly composed of gram eosin staining solution, methyl violet staining solution and gram iodine solution. It is a simple and cheap cellulose staining solution. After dyeing, the cellulose is blue.

### Protocols (for reference only)

1. Conventionally make paraffin section and dewax to water.
2. Add Gram Eosin Staining Solution, dye for 10 min and wash quickly with distilled water to avoid destaining.
3. Add Methyl Violet Staining Solution, dye for 3-5 min, and wash quickly with distilled water..
4. Add Gram Iodine Solution and dye for 3-5 min.
5. Discard Gram Iodine Solution and absorb the excess dye with absorbent paper.
6. Mix the same amount of aniline and xylene as differentiation solution and differentiate for about 30s.
7. Replace with new xylene to clean the tissue. Seal with resinene.

### Result

Cellulose	Blue
Back	Red

### Note

1. After staining with Gram Eosin Staining Solution, make sure there is enough red on the tissue after washing.
2. When treating with aniline xylene differentiation solution, pay attention to shaking gently to make the tissue decolorization uniform.
3. For your safety and health, please wear experimental clothes and disposable gloves.