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Cat: G2071

Size: 50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	3×50mL	Storage
Reagent(A1):HID Solution A	50mL	RT, avoid light
Reagent(A2):HID Solution B	3mL	RT, avoid light
Before use, mix A1 with A2 as the ratio of 50:3 to form HID Solution. It is not suitable to prepare in advance.		

Introduction

The principle of High Iron Diamine/Alcian Blue(HID/AB) Mucin Stain Kit is that the diammonium salt in the high iron diamine salt combines with sulfuric acid mucus to form a complex and be colored. When the pH value is much lower than 2.5, the sulfate radical in the tissue is ionized with a negative charge, forming a salt bond with the cation in HID, making the tissues with sulfate radical (such as sulfuric acid mucus) dye, sulfuric acid mucus material, etc. form brownish purple to brownish black.

HID Working Solution(G2070-A) is mainly used to identify sulfated acid mucus and sialic acid mucus. Nitrogen acetylated sialic acid mucus is produced in small intestine epithelium, oxygen acetylated sialic acid mucus and sulfated acid mucus are produced in large intestine epithelium. This method can be used with AB-PAS staining to identify the type of intestinal metaplasia and the type of mucus in metastatic tumors.

Protocols(for reference only)

1. Dewax the paraffin section to distilled water.
2. Soak in the prepared HID Solution for 12-24 hours (20-25 °C). Rinse with distilled water for 5mins.
3. (Optional)Add Alcian Staining Solution and stain for 10-20mins. Slightly wash with distilled water.
4. (Optional)Re-dyeing with Nuclear Fast Red Staining Solution for 5-10mins. Rinse in distilled water for 1min.
5. Dehydrate in series of ethanol, transparent by xylene and seal with resinene.

Result

Sulfated acid mucus (such as sulfated mucin and salivary mucin)	Brownish Purple to Brownish Black
Carboxylated protein acidic mucus substances (such as proteoglycans and hyaluronic acid)	Blue
Nucleus	Red

Note

1. The fixative adopts 10% neutral formalin.
2. HID Solution is not suitable to be prepared in advance. After use, it can be stored at 4 °C and used for 1-2 times, but its specificity is not good.
3. HID Solution is harmful to human body. Please operate it carefully to avoid contact with human skin.
4. For your safety and health, please wear experimental clothes and disposable gloves.