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Cat: G2101

Size: 500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

The fixative can be divided into aldehyde fixative, mercury fixative, alcohol fixative, oxidant fixative, picric acid salt fixative, etc. Formaldehyde in aldehydes and alcohol in alcohols are more commonly used.

Alcohol-Formalin Fixative is mainly composed of alcohol and formalin, which belongs to dehydration-crosslinking fixative. It can produce excellent results in the immunohistochemical reaction of identifying specific antibodies, and it is also conducive to the fixation of large samples with a large amount of fat.

Protocol(for reference only)

Operate according to the specific requirements of the experiment. Generally, the fixing time of specimens is 1-4h / mm, most of them are 4-24h. The fixing time of large specimens should be extended properly.

Note

1. Alcohol-Formalin Fixative has certain irritation and corrosiveness. Please operate carefully in ventilated environment
2. The fixing time changes due to different thickness of tissue samples. The suitable thickness of conventional biopsy tissue is 2-4mm, generally no more than 6mm. Proper material selection for the tissue is beneficial to the permeation of the fixative.
3. The effect of temperature on fixation is obvious. Increasing the temperature can accelerate the fixation, but the temperature should not be too high.
4. After fresh tissue is taken out, it should be fixed in time. If it can not be fixed in time, it should be stored in physiological saline for inspection in time.
5. For your safety and health, please wear experimental clothes and disposable gloves.