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Sperm Morphology Stain Kit(Shorr Method)

Cat:G2570

Size:3×20mL/3×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×20mL	3×100mL	Storage
Reagent A:Scarlet Orange Solution	20mL	100mL	RT, avoid light
Reagent B:Mordant Solution	20mL	100mL	RT, avoid light
Reagent C:Fast Green Solution	20mL	100mL	RT, avoid light

Introduction

The evaluation methods of sperm morphology recommended by WHO are Papanicolaou Method, Shorr Method and Diff-Quick Method.

The principle of Sperm Morphology Stain Kit(Shorr Method) is similar to that of Papanicolaou Method, because the proteins with different isoelectric points in sperm and cells have different charges under the same acidity, they can selectively combine with the corresponding dyes to dye. The nucleus is composed of acid substance, which has a strong affinity with basic dyes, while the cytoplasm has a strong affinity with acid dyes. It is especially suitable for the staining of sperm and cell samples such as pleural fluid, ascites and sputum.

Protocol(for reference only)

1. Fix the cell smear with 95% ethanol-glacial acetic acid (3:1) fixative for 10-15 min.
2. Soak in 80%, 70% and 50% ethanol for 1 min respectively.
3. Soak or wash in distilled water or tap water for 1min.
4. Stain with hematoxylin staining solution for 3-5min. Rinse with tap water for 2min.
5. (Optional) Differentiate by 0.5% hydrochloric acid ethanol differentiation solution for 4-5s. Rinse with tap water for 2min.
6. Blue in the bluing solution for 4 min.Rinse with tap water for 2min.
7. Stain with drops of Scarlet Orange Solution stain for 1-2 min.
8. Distilled water washed quickly, then treated with dropwise Mordant Solution for 1-2 min.
9. Drops of Fast Green Solution were stained for 4 min.
10. After pouring off the excess staining solution, 1% glacial acetic acid was added directly dropwise for 30s for differentiation.
11. Wash with water, and then dehydrate in 50%, 70%, 80%, 95% and absolute ethanol.
12. Transparent by xylene and seal with resinene.

Result

Acrosome Region of Sperm Head	Light Blue
Postacrosomal Region	Deep Blue
Middle Part	Slight Red
Rear Part	Blue or Slight Red
The excessive residual cytoplasm that usually located in the lower part of the head or around the middle part is dyed orange red	

Note

1. All dye solutions should be filtered and changed frequently.
2. For your safety and health, please wear experimental clothes and disposable gloves.