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Cat: G2581

Size: 10mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	10mL	Storage
Reagent A: Eosin Staining Solution	4mL	RT, avoid light
Reagent B: Nigrosin Staining Solution	6mL	RT, avoid light

Introduction

Sperm Vitality Kit (Eosin-Nigrosin Method) can be used to detect the survival rate of all sperm samples, and to identify the sperm with intact cell membrane by dye rejection method, so as to obtain the percentage of living sperm. The dye solution is composed of Eosin Stain Solution and Nigrosin Stain Solution. The dye rejection principle is that the damaged cell membrane, such as the membrane on the dead cell, allows the non membrane penetrating dye to enter for dyeing, while the cell membrane of the living cell can resist the dye entry, resulting in the phenomenon of dye rejection. The staining solution is only suitable for clinical diagnosis or other purposes.

Protocol (for reference only)

1. Put the fresh semen sample in 37 °C incubator or water bath for 30 min to liquefy the sperm, shake and mix well.
2. Take a clean small test tube, add 1 drop of fresh semen and 1-2 drops of Eosin Staining Solution, mix well, and place for 15s.
3. Add 3 drops of Nigrosin Staining Solution and mix well to form Semen-Eosin-Nigrosin Staining Solution, place for 30-60 s.
4. Add 1 drop of the above Semen-Eosin-Nigrosin Staining Solution on the slide to make a smear.
5. Dry it and view under the microscope. Under the oil immersion lens, the colorless sperm is white, the dead sperm is red and the background is reddish brown.
6. Count 200 sperm and calculate the percentage of unstained (living sperm) in 200 sperm.

Result

Living sperm	Colorless
Dead sperm	Red
Background	Reddish Brown

Note

1. Once the semen sample is liquefied, the sperm survival rate should be detected immediately, preferably within 30 min and no more than 1 h under any circumstances, so as to avoid causing sperm lose effect due to dehydration or temperature change and showing inaccurate staining results.
2. After the smear is dried, a small amount of salt crystals may be observed, which is a normal phenomenon.