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Sperm Vitality Detection Solution (Eosin Method)

Cat:G2582

Size:10mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Eosin is one of xanthene dyes in synthetic dyes, which is pink or pink powder. The molecular formula is $C_{20}H_6O_5Br_4Na_2$, and the molecular weight is 691.86. In clinical, the sperm survival rate can be evaluated by detecting the integrity of sperm membrane. The Sperm Vitality Detection Solution (Eosin Method) can be used to detect the survival rate of all sperm samples, and to identify the sperm with intact cell membrane by dye rejection method, so as to obtain the percentage of living sperm. The dye solution is composed of eosin and distilled water. The dye rejection principle is that the damaged cell membrane, such as the membrane on the dead cell, allows the non membrane penetrating dye to enter for dyeing, while the cell membrane of the living cell can resist the dye entry, resulting in the phenomenon of dye rejection. The staining solution is only suitable for clinical diagnosis or other purposes.

Protocol(for reference only)

1. Take a clean slide, add 1 drop of fresh semen and 1 drop of Sperm Vitality Detection Solution (Eosin Method), and mix well.
2. Cover and place it for 30-60s, then immediately view under high power microscope (such as 400 ×).
3. Count 200 sperm and calculate the percentage of unstained (living sperm) in 200 sperm.

Result

Living sperm	Colorless
Dead sperm	Red

Note

1. Once the semen sample is liquefied, the sperm survival rate should be detected immediately, preferably within 30 minutes and no more than 1 h under any circumstances, so as to avoid causing sperm lose effect due to dehydration or temperature change and showing inaccurate staining results.
2. For your safety and health, please wear experimental clothes and disposable gloves.