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Cat: G2590

Size: 50mL/100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit components

Reagent	50mL	100mL	Storage
Reagent A: S-M Solution A	1.5mL	3mL	2-8°C, avoid light
Reagent B: S-M Solution B	48.5mL	97mL	2-8°C, avoid light
Before use, mix Reagent A with B as the ratio of 3:97 to prepare S-M Working Solution. Use after filtering.			

Introduction

There are many methods to detect urinary sediment, such as quantitative examination, centrifugal precipitation examination, staining examination and so on. Urate Stain Kit(Sternheimer-Malbin Method) is a kind of staining urine sediment examination method. It is composed of crystal violet and safranin. After staining, the red blood cells are light purple and the white blood cells are orange red. The reagent is only used in scientific research and is not suitable for clinical diagnosis or other purposes.

Self Provided Materials

Incubator, slide, blood cell counting plate, cover glass, optical microscope

Protocols(for reference only)

1. Take out the reagent and restore to room temperature (25°C) for 30min.
2. Prepare S-M Working Solution and filter it for standby.
3. Take 0.2ml sediment, add 1 drop of S-M Working Solution and mix.
4. Cover the cover glass or fill the counting plate, place it for 3min and conduct microscopic examination.

Result

Red blood cell	Light Purple
Polymorphonuclear leukocyte	The nucleus is orange purple, and granules can be seen in the cytoplasm
Cell tube type	Deep Purple
Transparent tube type	Pink Red or Light Purple
Epithelial cells	The nucleus is purple and the cytoplasm is light purple or peach color

Note

1. Prepare a clean and dry urine collection cup. Generally reserve the middle urine.
2. In the process of centrifugation to obtain cell precipitation, for special cells, if the cell precipitation is insufficient, the centrifugation force can be appropriately increased or the centrifugation time can be prolonged.
3. In order to maintain the original morphological characteristics of cell components in urinary sediment, it is required to send them for examination quickly and complete the examination within 2h.
4. For your safety and health, please wear experimental clothes and disposable gloves.