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**Cat:** G2591

**Size:** 100mL

**Storage:** 2-8°C, avoid light, valid for 3 months.

### Introduction

There are many methods to detect urinary sediment, such as quantitative examination, centrifugal precipitation examination, staining examination and so on. Urate Stain Solution(Sternheimer Method) is a kind of staining urinary sediment examination method. It is composed of alician blue and pyronin . After staining, the red blood cells are red or colorless, and the white blood cells are dark blue, light blue or colorless. The reagent is only used in scientific research and is not suitable for clinical diagnosis or other purposes.

### Protocols(for reference only)

1. Take 0.2ml urinary sediment and add 1 drop(20 $\mu$ L) of Urate Stain Solution(Sternheimer Method) to mix well.
2. Cover the cover glass or fill the counting plate, place it and conduct microscopic examination.

### Result

Red blood cells	Red or Colorless
Polymorphonuclear leukocyte	Dark Blue, Light Blue or Colorless
Squamous epithelial cell	Light Pink or Purplish Red
Transitional epithelial cells, renal tubular epithelial cells	Purplish Red
Cell tube type	Light Blue or Dark Blue
Tubular particle	Pink or Purplish Red

### Note

1. Prepare a clean and dry urine collection cup. Generally, it is best to reserve the middle urine.
2. In the process of centrifugation to obtain cell precipitation, for special cells, if the cell precipitation is insufficient, the centrifugation force can be appropriately increased or the centrifugation time can be prolonged.
3. In order to maintain the original morphological characteristics of cell components in urinary sediment, it is required to send them for examination quickly and complete the examination within 2 hours.