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Electrolytic Decalcifying Solution(Formic Acid Method)

Cat: G2672

Size: 500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

In the process of tissue sectioning, when some tissues contain bone or calcification, the tissue containing calcium should not be directly embedded in paraffin. This is because the density between calcium and paraffin is different, it is difficult to cut a complete section. It is better to fix the calcium containing tissue before decalcification or conduct both at the same time. Then continue operations such as dehydration, transparency, wax immersion, embedding and slicing. There are many decalcification reagents, including organic acid, inorganic acid, EDTA and electrolytic decalcification. Among organic acids, formic acid and acetic acid are more suitable for decalcification of bone marrow.

The advantages of electrolytic decalcification solution are: ① rapid decalcification; ② suitable for most tissues. The disadvantages are: ① damage the tissue to some extent; ② it is need to carry out alkali treatment with sodium sulfate solution after decalcification.

Protocol(for reference only)

1. When the bone tissue is decalcified, pick up the material avoiding too thick, generally about 5mm.
2. After fixing the tissue, wash with PBS for three times.
3. Wash the tissue with distilled water for three times.
4. Installation of electrolysis device: take one platinum wire, wrap the sample with one end of platinum wire for several turns, place it on one side of the bottom of the container, and connect the other end of the platinum wire with the positive pole of the power supply. On the other side of the container, insert the lower end of a steel sheet or carbon rod into the container, and connect the upper section with the negative pole of the power supply.
5. Add Electrolytic Decalcification Solution and place it in the constant temperature water bath. Generally, the DC power supply is 6V, and if there is a rectifier, can adjust to 8-10V.
6. Rinse with running water for several hours.
7. Add 5% sodium sulfate solution for alkali treatment.
8. Conventional dehydration and embedding.

Note

1. The decalcification time of 5 mm thick bone tissue block is generally 10min-30 hours.
2. Decalcification should be thorough to prevent insufficient or excessive decalcification. The degree of decalcification should be controlled to shorten the decalcification time as much as possible without affecting the tissue section, so as to avoid tissue damage caused by too long decalcification.
3. It is better to fix the calcium containing tissue before decalcification or conduct both at the same time in order to reduce the degree of tissue damage.
4. The degree of decalcification should be detected in a while. Excessive decalcification will increase the degree of tissue damage and affect the staining results.
5. For your safety and health, please wear experimental clothes and disposable gloves.

Appendix:

Determination of the end point of decalcification (physical method): acupuncture, hand pinching, clamp and other methods are used to stop decalcification when the bone tissue becomes soft or there is no sense of resistance during acupuncture. Physical detection will damage the tissue structure to some extent, and try to avoid excessive force or repeated detection.